

# J4<sup>E</sup>

Group 2: Single Power (K0835)



## Features

- Standard 5 mph motors offer advanced performance, precise control and responsiveness via innovative single-stage gearbox
- Bold aesthetics
- Available with up to 10 inches of seat elevation

## Seating Options

- Q4 50° Tilt
- Power articulating foot platform
- 10 inch seat elevation

## Specifications

<b>DRIVE WHEELS</b>	14" Solid
<b>CASTER WHEELS</b>	<b>Front:</b> 6" Solid <b>Rear:</b> 6" Solid
<b>ANTI-TIP WHEELS</b>	N/A
<b>MAXIMUM SPEED</b> <sup>1</sup>	Up to 5 mph
<b>GROUND CLEARANCE</b>	3" (battery tray)
<b>TURNING RADIUS</b> <sup>2</sup>	20.5"
<b>OVERALL LENGTH</b> <sup>2</sup>	35.5"
<b>BASE WIDTH</b>	24.4"
<b>SEATING SIZES</b> (Q4 LIFT TILT)	<b>Width:</b> 16", 18" or 20" <b>Depth:</b> 16-22" <i>Coming soon: 22" width</i>
<b>SEAT-TO-FLOOR HEIGHTS</b> <sup>5</sup> (Q4 LIFT TILT)	17.5-19.5"
<b>BATTERY SIZE</b> <sup>3</sup>	NF-22
<b>BATTERY WEIGHT</b> <sup>7</sup>	38 lbs.
<b>AVAILABLE ELECTRONICS</b>	120A Q-Logic 3e <b>75A Q-Logic NE+ (Coming Soon)</b>
<b>BATTERY CHARGER</b>	8A Off-board
<b>MOTOR PACKAGES</b>	2-Pole Mid-Wheel 6 <sup>®</sup>
<b>WEIGHT CAPACITY</b> <sup>4</sup>	300 lbs.
<b>BASE WEIGHT</b>	138.6 lbs.
<b>MEDICARE CODE COVERAGE</b>	<b>Group 2: Single Power (K0835)</b>

## Available Colors



Red Corvette



Electric Blue



Back in Black

- 1) Speed varies with user weight, terrain type, battery charge/condition, and tire pressure.
- 2) Without front riggings
- 3) Battery capacity may change when adding power positioning systems.
- 4) Weight capacity is dictated by the selection of seating.
- 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 7) Battery weight may vary +/- 3 lbs. based on manufacturer.



quantumrehab.com

866.800.2002

**QUANTUM**<sup>®</sup>  
#1 FOR REHAB POWER