

**Jazzy® 1420 Series  
Group 3**  
Order Form / Price List  
Revised: February 2011



**Pride Mobility Products Corporation®**  
182 Susquehanna Ave., Exeter, PA 18643  
Sales: Phone: (800) 800-8586 • Fax: (800) 800-1636

This product is available for purchase online at [www.mypride.com](http://www.mypride.com)

Account #:	Date:	Order No.	Fax No.	Phone No.
Bill To:	Ship To:		P.O. #:	
Address:	Address:		Discount:	
Phone Contact:	Marked For:		Terms:	
			Freight:	
			Taken By:	

DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement. HCPCS codes should not be considered as legal advice.

All pricing is MSRP.

**BASE MODEL**

- Jazzy 1420 3VHD-C K0853 \$8,700
- Jazzy 1420 3VHD-SS K0852 \$8,700

**COLOR**

- Champagne NC
- Candy Apple Red NC
- Viper Blue NC
- Onyx Black NC
- Sunburst Orange NC

**BATTERIES (Installed)**

- Group 24, 70 amp, Gel E2363 \$352ea

**BATTERY CHARGER**

- 8 amp Off-board Charger Std.

**NON-EXPANDABLE CONTROLS**

- PG VR2 70 amp<sup>1</sup> Std
- PG 6-key VR2 70 amp<sup>1,4</sup> \$200

**EXPANDABLE CONTROLS**

All joysticks must be used in conjunction w/ a PG controller

- CONTROLLER OPTIONS**
- PG Controller \$695

All joysticks or enhanced display must be used in conjunction w/ a PG controller

- Harness for Expandable Controller \$350

- JOYSTICK OPTIONS**
- PG Electronics**
- 3-Key Remote Plus 70 amp<sup>1,3</sup> \$700  
(Select Box for Profile)
    - 1 Profile  Multiple Profile

- Remote Plus 70 amp with Speed Dial<sup>3,5</sup> \$700

- Remote Plus 70 amp with Toggles<sup>3,5</sup> \$700

- Remote Plus 70 amp with Switch Jacks<sup>3,5</sup> \$700

- Microdrive, 70 amp<sup>3,5,6</sup> \$700

- Omni Plus<sup>17</sup> \$1721.63  
(Requires a secondary control to operate the power chair)

Joystick includes electronics & junction boxes

Pride has submitted for the establishment of a new HCPCS code for the Upgraded Proportional Joystick. Call for clarification.

**PROGRAMMER**

- PG Programmer \$423

**JOYSTICK MOUNT**

- Fixed, In-Line Joystick Mount Std.
  - Left  Right
- Swing-away Joystick Mount<sup>18</sup> E1028 \$210
  - 4" <sup>19</sup>  7"
  - Left  Right
- Multi-axis Joystick Mount K0108 \$100
  - Left  Right  
(Fixed, does not swing-away)
- Vertical Adj. Clamp<sup>8</sup> K0108 \$95  
Provides vertical joystick adjustment w/o adjusting armrest height.

**JOYSTICK HANDLE OPTIONS**

- PG Controls**
- Std. Knob Std
  - Mushroom Handle E2323 \$71
  - T-Bar E2323 \$107
  - Chin Cup E2324<sup>15</sup> \$71
  - Extended I-Handle E2323 \$71
  - Soft Ball Handle E2323 \$50
  - Medium Ball Handle E2323 \$50
  - Goal Post Handle E2323 \$66
    - 3" W  4" W

# Jazzy® 1420 Series Group 3

Order Form / Price List  
Revised: June 2011

## DRIVE WHEELS

<input type="checkbox"/> 16" Pneumatic (pair)	Std
<input type="checkbox"/> Flat-Free Inserts (pair)	\$213

## CASTER WHEELS

<input type="checkbox"/> Pneumatic (pair)	Std
<input type="checkbox"/> Flat-Free Inserts (pair)	\$105

## NEW SWING-AWAY LEGRESTS

<input type="checkbox"/> 70o Swing-away footrests Style #8	\$162
<input type="checkbox"/> 80o Swing-away footrests Style #8	\$162
<input type="checkbox"/> 70o Tapered S/A footrests Style #8	\$162

(N/C with power positioning systems)

<input type="checkbox"/> Quad Release Lever	NC
---	----

<b><input type="checkbox"/> S/A Standard Extensions:</b>	NC
--	----

<input type="checkbox"/> 11"-13" Ext.	NC
<input type="checkbox"/> 13"-16.5" Ext	Std.
<input type="checkbox"/> 15.5"-19.5" Ext.	NC
<input type="checkbox"/> 18.5"-22" Ext.	NC

(Lengths over 17" may cause interference)

<input type="checkbox"/> Extruded Footplate	Std.
<input type="checkbox"/> 7.25" x 6" 10 (Tapered legs require 18"W or greater seat)	
<input type="checkbox"/> 8.5 x 6" 9 (Tapered legs require 20"W or greater seat)	

### Angle Adjustable Footplates (pair) K0040 \$150

<input type="checkbox"/> 4"x6" 16 (Tapered legs require 16"W or greater seat)	
<input type="checkbox"/> 5"x6" 10 (Tapered legs require 18"W or greater seat)	
<input type="checkbox"/> 5"x8" 10 (Tapered legs require 18"W or greater seat)	
<input type="checkbox"/> 6"x8" 9 (Tapered legs require 20"W or greater seat)	
<input type="checkbox"/> 6"x11" 9 (Tapered legs require 20"W or greater seat)	
<input type="checkbox"/> 8"x11" 11 (Tapered legs require 24"W or greater seat)	
<input type="checkbox"/> Wheel Bumpers	\$25

### Multi-Axis Ext. (ball & socket):

<input type="checkbox"/> 11"-16.5" Ext.	
<input type="checkbox"/> 15.5"-19.5" Ext.	
<input type="checkbox"/> 18.5"-22" Ext.	

(Lengths over 17" may cause interference)

### Multi-Axis Footplates (pair) K0040 \$495

<input type="checkbox"/> 4"x6" 16 (Tapered legs require 16"W or greater seat)	
<input type="checkbox"/> 5"x6" 10 (Tapered legs require 18"W or greater seat)	
<input type="checkbox"/> 5"x8" 10 (Tapered legs require 18"W or greater seat)	
<input type="checkbox"/> 6"x8" 9 (Tapered legs require 20"W or greater seat)	
<input type="checkbox"/> 6"x11" 9 (Tapered legs require 20"W or greater seat)	
<input type="checkbox"/> 8"x11" 11 (Tapered legs require 24"W or greater seat)	
<input type="checkbox"/> Wheel Bumpers	\$25

## FRONT RIGGINGS

<b>Rigging Type (pair)</b>	<b>E0990</b>	<b>\$325</b>
<input type="checkbox"/> Elevating Leg Rests		
Select footplate size for above:		
<input type="checkbox"/> 71/4 x 6	<input type="checkbox"/> 81/2 x 6	
Adjustable Calf Pads w/ ELRS (1 hole std.)		
<input type="checkbox"/> 1 hole	<input type="checkbox"/> 3 holes*	<input type="checkbox"/> 6 holes*
<b>*Potential for interference w/ shrouds</b>		
<b>Extensions:</b>		
<input type="checkbox"/> 11" - 14" Ext.		
<input type="checkbox"/> 14" - 17" Ext.		
<input type="checkbox"/> 17" - 20" Ext.		
(Lengths over 17" may cause interference)		

<input type="checkbox"/> HD Drop-Ins, 70 <sup>o2,7,13</sup>	\$420
<input type="checkbox"/> HD Drop-Ins, 80 <sup>o2,7,13</sup>	\$420
<b>N/A w/ contoured or Specialty seating</b>	
<b>Angle Extensions:</b>	NC
<input type="checkbox"/> 12" - 14.5" Ext.	
<input type="checkbox"/> 14.5" - 18" Ext.	
<input type="checkbox"/> 17.5" - 21" Ext.	
(Lengths over 17" may cause interference)	

<input type="checkbox"/> HD Foot Platform (Substitute for S/A)	NC
<input type="checkbox"/> No Front Rigging	NC

### Footplates (pair)

<input type="checkbox"/> Reinforced Footplate	Std.
<input type="checkbox"/> 7.25" x 6" <sup>10</sup>	<input type="checkbox"/> 8.5 x 6" <sup>9</sup>
<input type="checkbox"/> 5"x8" Angle-Adj. Footplate <sup>10</sup>	<b>K0040</b> \$150
<input type="checkbox"/> 6"x8" Angle-Adj. Footplate <sup>9</sup>	<b>K0040</b> \$150
<input type="checkbox"/> 6"x11" Angle-Adj. Footplate <sup>9</sup>	<b>K0040</b> \$150
<input type="checkbox"/> 8"x11" Angle-Adj. Footplate <sup>11</sup>	<b>K0040</b> \$150
<input type="checkbox"/> Wheel Bumpers (Ang-adj. only)	\$25

### Extensions:<sup>12</sup>

<input type="checkbox"/> 11"-14" Ext.	NC
<input type="checkbox"/> 14"-17" Ext	Std.
<input type="checkbox"/> 17"-20" Ext.	NC

<input type="checkbox"/> High mount Clamp-on; Ang. Adj.	<b>K0040</b> \$326
<input type="checkbox"/> Small	<input type="checkbox"/> Large
<input type="checkbox"/> High mount Clamp-on; Multi-axis	<b>K0040</b> \$412
<input type="checkbox"/> Small	<input type="checkbox"/> Large

<input type="checkbox"/> High mount footplatform <sup>19</sup>	<b>K0108</b> \$450
<input type="checkbox"/> HD Foot Platform	
Lower Extensions	
<input type="checkbox"/> 10"-14" L	
<input type="checkbox"/> 14"-17" L	
<input type="checkbox"/> 17"-20" L	

### Rigging Options

<input type="checkbox"/> Heel Loops (pair) <sup>13</sup>	<b>E0951</b> \$38
<input type="checkbox"/> Gel Neoprene Footrest	<b>K0108</b>
Sleeve (pair) <sup>14</sup>	\$65
<input type="checkbox"/> ELR/ALR Adductor Buttons (pair)	<b>K0108</b> \$100
<input type="checkbox"/> POSAlinc Knee Adductor (pair)	
for S/A legrests	<b>K0108</b> \$464
<input type="checkbox"/> Residual limb support left	<b>E1020</b> \$300
<input type="checkbox"/> Residual limb support right	<b>E1020</b> \$300

## Footnotes for pages 1 & 2

- Std. w/ one programmable drive profile.
- May cause interference w/ front casters w/ certain seat widths and lower leg lengths. Call Customer Service.
- Needed for controlling power positioning systems.
- For use with power elevating seat.
- Avail. w/ five individual drive profiles.
- Requires multi-axis joystick mount; N/C for mount.
- N/A w/ contoured seating.
- N/A on Specialty, contoured seating or waterfall armpads.
- Only avail. on 16" W or greater seats.
- Only avail. on 14" W or greater seats.
- Only avail. on 20" W or greater seats.
- Knee to heel length is measured from the top of the legrest to the footplate. It will increase when contoured or rehab seating w/ cushions is selected. The overall length may vary approx. 3" depending on choice of seat and/or cushion.
- Std. w/ ang-adj. footplates.
- Avail. w/ S/A footrests & 70 degree HD drop-ins.
- HCPCS code E2323 when not used in conjunction with chin control.
- Only avail. on 12" W or greater seats.
- Please consult specialty controls order form for a variety of secondary controls and mounting options.
- Can only mount one armrest accessory in combination w/ S/A joystick.
- N/A with waterfall armpads and dual post armrests.

# Jazzy® 1420 Group 3

Order Form / Price List  
Revised: June 2011

## SEATING

High-Back w/ Headrest Std.

- 22" W x 20" D
- 24" W x 20" D
- 24" W x 22" D

Gray Vinyl  Beige Vinyl

High-Back w/ Headrest  
& Solid Seat Pan<sup>1</sup> Std.

- 22" W x 20" D
- 24" W x 20" D
- 24" W x 22" D

Gray Vinyl  Beige Vinyl

**22" W AND ABOVE WILL NOT BE LTD RECLINE.**

Custom seat dimensions are available.  
Please contact the Custom Group for pricing and leadtime.

## ARMRESTS

For High-Back Seats

- Flip-up, removable Std.
- Flip-up, height-adjustable \$170  
removable<sup>4</sup> E0973

Armpads

- Desk Length, 10" NC
- Left NC  Right NC
- Full Length, 14" NC
- Left Std.  Right Std.

## POWER ELEVATING SEAT<sup>4,6,7</sup>

(Max weight capacity 300lbs.)

- Power Elevating Seat \$975
- Through 6-Key VR2 NC
- Power Elevating Seat \$1375  
Through Remote Controller
- Lap Belt Std.

## MANUAL RECLINER<sup>2,5</sup>

Manual Recline High Back E1226 \$600

Manual Recline High Back<sup>1</sup>  
w/ Solid Seat Pan E1226 \$600  
(180 degrees of Recline)  
(Dimensions 18"W x 18"D)

- |  |  |
|--|--|
| <input type="checkbox"/> Gray Vinyl                                  | <input type="checkbox"/> Beige Vinyl                                 |
| <b>Width:</b>  | <b>Depth:</b>  |
| <input type="checkbox"/> 16" <span style="float: right;">NC</span>   | <input type="checkbox"/> 16" <span style="float: right;">NC</span>   |
| <input type="checkbox"/> 18" <span style="float: right;">Std.</span> | <input type="checkbox"/> 18" <span style="float: right;">Std.</span> |
| <input type="checkbox"/> 20" <span style="float: right;">NC</span>   | <input type="checkbox"/> 20" <span style="float: right;">NC</span>   |

When ordering a lap belt with this seat, an additional 10" is added to accommodate usable length

## CUSHIONS<sup>3</sup>

- No Cushion Required
- Simplicity Cushion E2601  
16 x 16 - 20 x 20 \$98
- Solution 1 Cushion
- 14 x 14 - 20 x 18 E2603 \$380
- 20 x 20 - 20 x 24 E2603 \$430
- 22 x 16 - 24 x 24 E2604 \$430
- Solution Cushion
- 10 x 10 - 20 x 18 E2607 \$380
- 20 x 20 - 20 x 24 E2607 \$430
- 22 x 16 - 24 x 24 E2608 \$430
- Spectrum Cushion
- Gel
- 10 x 10 - 20 x 18 E2605 \$480
- 20 x 20 - 20 x 24 E2605 \$530
- 22 x 16 - 24 x 24 E2606 \$530
- Visco
- 10 x 10 - 20 x 18 E2607 \$480
- 20 x 20 - 20 x 24 E2607 \$530
- 22 x 16 - 24 x 24 E2608 \$530
- Spectrum Air
- 14 x 12 - 20 x 24 E2622 \$480
- 22 x 16 - 24 x 24 E2623 \$530
- Spectrum Air Contour
- 14 x 12 - 20 x 24 E2624 \$499
- 22 x 16 - 24 x 24 E2625 \$549

### Footnotes for pages 3:

1. Includes black vinyl seat pan, back available in all listed colors.
2. 18" x 18" std. Special measurements are not available.
3. See Synergy Cushion order forms for detailed pricing, sizes & codes.
4. N/A on manual recliner.
5. Rear mount accessories are N/A.
6. Includes appropriate actuator module.
7. Max. weight capacity is 300 lbs.

# Jazzy® 1420 Group 3

Order Form / Price List  
Revised: June 2011

## Synergy Seating Systems

Synergy Seating™ Frame **\$434**

### WIDTH

#### Set Seat Width At:

<input type="checkbox"/> 10"	NC	<input type="checkbox"/> 11"	NC
<input type="checkbox"/> 12"	NC	<input type="checkbox"/> 13"	NC
<input type="checkbox"/> 14"	NC	<input type="checkbox"/> 15"	NC
<input type="checkbox"/> 16"	NC	<input type="checkbox"/> 17"	NC
<input type="checkbox"/> 18"	Std.	<input type="checkbox"/> 19"	NC
<input type="checkbox"/> 20"	NC	<input type="checkbox"/> 21"	NC
<input type="checkbox"/> 22"	NC	<input type="checkbox"/> 23"	NC
<input type="checkbox"/> 24"	NC		

#### Seat Width Package

<input type="checkbox"/> 10"W-14"W	<input type="checkbox"/> 14"W-18"W
<input type="checkbox"/> 16"W-20"W	<input type="checkbox"/> 20"W-24"W

### DEPTH

#### Set Seat Depth At:

<input type="checkbox"/> 10"	NC	<input type="checkbox"/> 11"	NC
<input type="checkbox"/> 12"	NC	<input type="checkbox"/> 13"	NC
<input type="checkbox"/> 14"	NC	<input type="checkbox"/> 15"	NC
<input type="checkbox"/> 16"	Std.	<input type="checkbox"/> 17"	NC
<input type="checkbox"/> 18"	NC	<input type="checkbox"/> 19"	NC
<input type="checkbox"/> 20"	NC	<input type="checkbox"/> 21"	NC
<input type="checkbox"/> 22"	NC	<input type="checkbox"/> 23"	NC
<input type="checkbox"/> 24"	NC		

#### Seat Depth Package

<input type="checkbox"/> 10"D-14"D	<input type="checkbox"/> 14"D-18"D
<input type="checkbox"/> 16"D-20"D	<input type="checkbox"/> 20"W-24"W

Please consult your Quantum Specialist or contact the Custom Group before ordering

### BACK CANES<sup>14</sup>

<input type="checkbox"/> Back Canes w/ push handles	Std.
<input type="checkbox"/> Back Canes w/out push handles	\$125
<input type="checkbox"/> Offset (Meso) \$195	<input type="checkbox"/> Inset (Endo) \$195
<input type="checkbox"/> 2" (1" per side)	<input type="checkbox"/> 4" (2" per side)
N/A with TRU-Comfort, Dual Post, Flip-back armrests may interfere with offset armrests .	
N/A w/ any reclining systems.	

Actual Back Cane Height is reflected in measurement

#### Back Cane Height

<input type="checkbox"/> 14"	NC	<input type="checkbox"/> 15"	NC
<input type="checkbox"/> 16"	NC	<input type="checkbox"/> 17"	NC
<input type="checkbox"/> 18"	Std.	<input type="checkbox"/> 19"	NC
<input type="checkbox"/> 20"	NC	<input type="checkbox"/> 21"	NC
<input type="checkbox"/> 22"	NC	<input type="checkbox"/> 23"	NC
<input type="checkbox"/> 24"	NC	<input type="checkbox"/> 25"	NC

### BACK OPTIONS

<input type="checkbox"/> Synergy Back up to 20" W	\$275
<input type="checkbox"/> Synergy Back 21-24" W	\$325
<input type="checkbox"/> Omit Back	
<b>Back Height:<sup>1</sup></b>	
<input type="checkbox"/> 18"H	<input type="checkbox"/> 19"H
<input type="checkbox"/> 20"H	<input type="checkbox"/> 21"H
<input type="checkbox"/> 22"H	<input type="checkbox"/> 23"H
<input type="checkbox"/> 24"H	<input type="checkbox"/> 25"H
<input type="checkbox"/> Sling Back	\$98

Back height will be the same as back cane height. Includes gap for cushion.

Push Bar selection must match back width; center tube standard when selecting flip-up arms.

### BACK OPTIONS CONT.

<input type="checkbox"/> TRU-Comfort Back <sup>9,13</sup>	
(Even sizes only 16" W - 22"W)	
<input type="checkbox"/> Std.	<input type="checkbox"/> New Plus
<b>Back Height:</b>	<input type="checkbox"/> 22" <input type="checkbox"/> 25"
<input type="checkbox"/> Back up to 20" W	<b>E2620</b> \$595
<input type="checkbox"/> 22" W	<b>E2621</b> \$695
<input type="checkbox"/> TRU-Comfort Cushion	
16 x 16 - 20 x 18	<b>E2601</b> \$380
20 x 20 - 20 x 22	<b>E2601</b> \$380
22 x 16 - 22 x 24	<b>E2602</b> \$480
<input type="checkbox"/> TRU-Comfort Plus Cushion	
16 x 16 - 20 x 18	<b>E2603</b> \$380
20 x 20 - 20 x 22	<b>E2603</b> \$380
22 x 16 - 22 x 24	<b>E2604</b> \$480
<input type="checkbox"/> Omit Back	-\$100
See TRU-Comfort Back & Cushion order forms for specific sizes	

TRU-Comfort Plus-utilizes a 4-way stretch material cover on cushions & backs

### ARMRESTS

<input type="checkbox"/> Flip-back, cane mounted <sup>2</sup>	Std.
<input type="checkbox"/> 2-Post, flip-back, height <sup>2</sup>	<b>E0973</b> \$150ea
adjustable, removable armrests (9"-13" of height adjustment in 1" increments)	
<input type="checkbox"/> Single Post height adj. <sup>7,8</sup>	\$75ea
quick height, removable <b>E0973</b>	
11"-17" of height adj. 3" on top in 1/2" increments & 3" on bottom in 1" increments.	
<input type="checkbox"/> Quad Release Lever	NC
<input type="checkbox"/> Removable, height adjustable	<b>E0973</b> \$75ea
heavy duty armrests <sup>12</sup>	
<input type="checkbox"/> 8"H-12"H	<input type="checkbox"/> 10"H-14"H
<input type="checkbox"/> Cantilever Armrests	SOF
Angle-Adjustable: - 30° through + 30°	
Height Adjustment 11" - 17" depends on back cane height. N/A w/ manual recline, flip-forward back canes, Cane Clamp Push Bar, Shape Back, Sling back & Specialty Seat Call for details when selecting aftermarket arm troughs or Universal Swivel w/ desk length arms.	

### ARMPADS

<input type="checkbox"/> Full length, Straight, 2" X 14"	Std.
<input type="checkbox"/> Left	<input type="checkbox"/> Right
<input type="checkbox"/> Desk Length, Straight, 2" X 10"	Std.
<input type="checkbox"/> Left	<input type="checkbox"/> Right
<input type="checkbox"/> Full Length, Waterfall, 2" X 14"	NC
<input type="checkbox"/> Left	<input type="checkbox"/> Right
<input type="checkbox"/> Desk Length, Waterfall, 2" X 10" 4,5	NC
<input type="checkbox"/> Left	<input type="checkbox"/> Right

### HEAD REST

<input type="checkbox"/> Stealth	<b>E0955</b>
<input type="checkbox"/> Large (11"W x 6" H)	<input type="checkbox"/> Small (9 3/4" W x 3" H) \$205
<input type="checkbox"/> Stealth Comfort Plus	<b>E0955</b>
<input type="checkbox"/> Small (10" W)	<input type="checkbox"/> Large (14" W) \$205
<input type="checkbox"/> Removable Head Rest Bracket	<b>E1028</b> \$210
<input type="checkbox"/> Stealth, swing-away (Flip-Down)	
Head Rest Bracket	<b>E1028</b> \$210
Interfaces w/ any headrest that uses 1/2" square down-tube.	
<input type="checkbox"/> Medical Applications	<b>E0955</b>
<input type="checkbox"/> Large (11"W x 6" H)	<input type="checkbox"/> Small (9 3/4" W x 3" H) \$230
<input type="checkbox"/> Removable Head Rest Bracket	<b>E1028</b> \$235
<input type="checkbox"/> Fixed, Offset Mounting Bracket	NC
<input type="checkbox"/> TRU-Comfort	<b>E0955</b>
<input type="checkbox"/> 8" W	<input type="checkbox"/> Leatherette \$205
<input type="checkbox"/> 10" W	<input type="checkbox"/> Plus \$205
<input type="checkbox"/> Cobra Removable Head Rest Bracket	<b>E1028</b> \$210
<input type="checkbox"/> Removable Bracket for	<b>E1028</b> \$210
after market headrest	
<input type="checkbox"/> No Headrest Required	

### ACCESSORIES

<input type="checkbox"/> Rear Basket <sup>10,11</sup>	\$103
<input type="checkbox"/> Cane/Crutch Holder <sup>10,11</sup>	<b>E2207</b> \$98
<input type="checkbox"/> Oxygen Holder <sup>10,11</sup>	<b>E2208</b> \$103
<input type="checkbox"/> Oxygen Holder <sup>10,11</sup>	<b>E2208</b> \$150
(Synergy Seat only)	
<input type="checkbox"/> Walker Holder <sup>10,11</sup>	<b>K0108</b> \$92
<input type="checkbox"/> Flag Holder (SOFO Power Positioning)	\$17
<input type="checkbox"/> Weather Cover, Large, Black	\$202
<input type="checkbox"/> Saddle bag(contoured seating only)	\$28
<input type="checkbox"/> Cup holder	\$17
<input type="checkbox"/> Left	<input type="checkbox"/> Right
N/A on desk length arm with S/A mount or toggle box.	
<input type="checkbox"/> Lap Belt	Std.
<input type="checkbox"/> 50" <sup>11</sup>	<input type="checkbox"/> 60" <input type="checkbox"/> 70" <input type="checkbox"/> 80" <input type="checkbox"/> 90"

### Footnotes for pages 4:

1. Measured from seat pan to top of back pan w/ 3" gap to accommodate the cushion.
2. S/A laterals interfere w/ flip-back arms.
3. Tools are required to obtain additional 4" of height adjustment at bottom of armrest. 3" of adjustment at top requires no tools. Total available adjustment range is 11" - 17".
4. Only avail. on removable, adjustable, quick-height arms.
5. N/A w/ S/A joystick mount.
6. Available w/ solid backs only.
7. Desk length armrests w/ an accessory limit the accessibility to the height adj. lever.
8. Tools are required to obtain additional 4" of height adjustment at bottom of armrest. 3" of adjustment at top requires no tools. Total available adjustment range is 11" - 17".
9. Cane mounted, flip-back arms are N/A.
10. Includes mounting bracket and hardware.
11. N/A on manual recliner high-back.
12. Desk length armrests w/ an accessory limit the accessibility to the height adj. lever.
13. Back height is measured from the top of a 3" cushion on the seat.
14. Height is measured from seat pan to top of cane.

#### Lowest seat-to-floor height

Synergy Seat:	18.125"
High-back Seat:	22.5"

All measurements +/- 1/8" with 18" x 18" seat, fixed towers (please see Tilt information above), 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Increased widths may increase seat-to-floor height. Contact customer service for details. Listed measurements are in the lowest position. It does not include elevating systems.

### Quantum Custom Delivery Service

Includes: liftgate service and delivery of your product(s) inside your facility Contact your Quantum Representative for details.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications and pricing without prior notice.