



## Quantum Outback Order Form with Contoured Seating

Account Number: \_\_\_\_\_ Date: \_\_\_\_\_  
 Provider Name: \_\_\_\_\_  
 Contact: \_\_\_\_\_  
 Phone: \_\_\_\_\_ Fax: \_\_\_\_\_  
 Email: \_\_\_\_\_  
 PO Number: \_\_\_\_\_  
 Quote Name/Marked for: \_\_\_\_\_  
 Ship to Address: \_\_\_\_\_  
 City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_



(TRU Balance Power Positioning shown above)

### INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.


This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email ([quantumorders@pridemobility.com](mailto:quantumorders@pridemobility.com)). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

**HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.**

Thank you for choosing Quantum!

### PATIENT INFORMATION

Weight: \_\_\_\_\_ Height: \_\_\_\_\_ A. Shoulder Width: \_\_\_\_\_ D. Max Sitting Height: \_\_\_\_\_ G. Thigh Depth: \_\_\_\_\_  
 B. Chest Width: \_\_\_\_\_ E. Shoulder Height: \_\_\_\_\_ H. Lower Leg Length: \_\_\_\_\_  
 C. Hip Width: \_\_\_\_\_ F. Axilla Height: \_\_\_\_\_ I. Elbow Height: \_\_\_\_\_



\_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_

### 1. BASE MODEL

**REQUIRED.** 350 lb. weight capacity with Q-Captain seating, 14" pneumatic drive wheels, base mounted lights, and 8 amp off board charger. Dimensions of base without legrests: 28.82"Wx42"L.

Quantum Outback ..... \$19,595  
 Part: Q OUTBACK. HCPCS: K0899 - Power mobility device, not coded by PDAC

**2. BASE COLOR**

**REQUIRED.** Select one shroud color.

Cadet Base Shroud Color ..... No Charge Part: DGN2301319-TGY3-PCL2	Desert Camouflage Base Shroud Color ..... <i>Coming soon</i> Part: DGN2301319-WCM11-PCL2
Tan Base Shroud Color ..... No Charge Part: DGN2301319-TTN1-PCL2	Tropic Lightning Camouflage Base Shroud Color ..... <i>Coming soon</i> Part: DGN2301319-WCM12-PCL2
Black Camouflage Base Shroud Color ..... <i>Coming soon</i> Part: DGN2301319-WCM9-PCL2	Jungle Camouflage Base Shroud Color ..... <i>Coming soon</i> Part: DGN2301319-WCM13-PCL2
White Camouflage Base Shroud Color ..... <i>Coming soon</i> Part: DGN2301319-WCM10-PCL2	

**3. MOTORS**

**REQUIRED.**

5 MPH Motors ..... No Charge Part: MOT2205971	6 MPH Motors ..... No Charge <i>Coming soon.</i> Part: MOT2300107
--	--

**4. POWER ADJUSTABLE SEAT HEIGHT**

**Optional.** ^Please verify the appropriate billing code to be used for this item with all third-party payors.

10" Power Lift for Captain Seat ..... \$3,725 Part: XRF000353. HCPCS: E2298/E2300^ 10" of scissor lift. 22"W x 22"D maximum seat size.
--

**5. OMIT MOUNTED LIGHTS**

**Optional.**

Omit Base Mounted Lights ..... No Charge Part: DGN2400803
--

**6. QUANTUM BACKUP CAMERA**

**Optional.** The backup camera may interfere with armrest mounted accessories such as media device holder.

Quantum Backup Camera ..... \$368 Part: MEC1910365
---

**1. Optional Camera Mounting Bracket**

Backup Camera Captain's Headrest Bracket ..... \$21 Part: MEC1901231	Backup Camera Track Mount Bracket ..... \$27 Part: MEC2104227
---	--

**2. Optional Display Mounting Bracket**

Backup Camera Q-Logic Joystick Mount Bracket ..... \$105 Part: MEC2106123
--



**7. BATTERIES**

**Optional.** The powerchair requires 2 batteries and will be installed if ordered.

Group 27, Interceptor 2, Gel Battery ..... \$519 ea. Part: BAT2401900
--

**8. JOYSTICK**

**REQUIRED.** Select a type of drive control from the following sections.

**8.1a. Omit Joystick**

No Joystick, Setup for Q-Logic Electronics Part: XRF000010 Requires selection of Enhanced Display
---

**8.1b. Q-Logic 3e Joystick**

Q-Logic 3e Joystick ..... No Charge Part: CTL1910600 Can operate full light package. Not available with alternative or specialty controls.
--

**8.1c. Expandable Controls**

Expandable controls allow multiple or specialized controllers.

**1. Q-Logic 3 Required Components**

Q-Logic 3 controls require the controller module and expandable harness below. If a Q-Logic joystick or specialty control is selected, Customer Service will quote the order with these items.

Q-Logic 3 EX Controller Module ..... \$730 Part: XRF000002. HCPCS: E2377	Q-Logic 3 Expandable Controller Harness ..... \$458 Part: XRF000264. HCPCS: E2313
---	--

**2. Q-Logic 3 Joystick**


Q-Logic 3 EX Joystick ..... \$840 Part: CTL167724	Q-Logic 3 EX Joystick for Full Light Package ..... \$840 Part: CTL167725 Required with full light packages.
--	---



**9. JOYSTICK HANDLES**

**Optional.** The standard Q-Logic 3 joystick knob can be seen in the image on the previous page. An optional handle selected below will be mounted to the joystick and the standard handle will be shipped with the unit.

Mushroom Handle (B) ..... \$90 Part: ACCASMB1860. HCPCS: E2323	Extended I-Handle (G) ..... \$90 Part: ACCASMB1857. HCPCS: E2323
Chin Cup (C) ..... \$90 Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control)	Bodypoint 3" Goal Post Handle ..... \$104 Part: PLSKNCP1057. HCPCS: E2323
T-Bar (D) ..... \$126 Part: ACCASMB1858. HCPCS: E2323	Bodypoint 4" Goal Post Handle ..... \$104 Part: PLSKNCP1058. HCPCS: E2323
Medium Ball Handle (E) ..... \$69 Part: ACCASMB1859. HCPCS: E2323	Q-Logic 3 Carrot Handle (A) ..... No Charge Part: DGN124308
Soft Ball Handle (F) ..... \$80 Part: ACCASMB1856. HCPCS: E2323	




**10. JOYSTICK MOUNTING BRACKETS**

**Optional.** Select a mount to be used with a joystick selected in previous section.

**10.1a. Fixed, In-line Joystick Bracket**


Fixed, In-line, Bracket ..... Standard Part: FRMASMB10617
--



**10.1b. Swing-Away (Retractable) Joystick Bracket**

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

Swing-away, Q-Logic 3 ..... \$319 Part: ACC1908434. HCPCS: E1028
---



**10.1c. Basic Swing Away Joystick Mount**

Basic Swing Away, Q-Logic, Left Side ..... \$245 Part: ACC2201041. HCPCS: E1028	Basic Swing Away, Q-Logic, Right Side ..... \$245 Part: ACC2201042. HCPCS: E1028
--	---

**10.1d. Specialty Joystick Mounts**

Select a bracket and side to mount it.

Retract4 Mount ..... \$368 Part: XRF000003. HCPCS: E1028 Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.	Multi-Axis Mount ..... \$126 Part: XRF000004 Allows height adjustment of joystick.
--	--



**10.2. Joystick Mount Location**

Left Side Mount Part: PTOINDV1225	Right Side Mount Part: PTOINDV1226
--------------------------------------	---------------------------------------

**11. CONTOURED SEATING**

**REQUIRED.**

**11.1. Occupied Transit, WC-19**

WC-19 Occupied Transit Option for Contoured Seat ..... \$368 <b>Coming soon.</b> Part: KIT122609 To determine if your chair is WC-19 compliant (occupied), please refer to this website: <a href="https://wc-transportation-safety.umtri.umich.edu/">https://wc-transportation-safety.umtri.umich.edu/</a>
--

**11.2a. Quantum Captain Seat**

16"W x 16"D ..... No Charge Part: SET2005145	18"W x 20"D ..... No Charge Part: SET2005151	20"W x 22"D ..... No Charge Part: SET2005909
16"W x 18"D ..... No Charge Part: SET2005146	18"W x 22"D ..... No Charge Part: SET2005908	22"W x 18"D ..... No Charge Part: SET2005155
16"W x 20"D ..... No Charge Part: SET2005147	20"W x 16"D ..... No Charge Part: SET2005152	22"W x 20"D ..... No Charge Part: SET2005156
18"W x 16"D ..... No Charge Part: SET2005149	20"W x 18"D ..... No Charge Part: SET2005153	22"W x 22"D ..... No Charge Part: SET2005157
18"W x 18"D ..... No Charge Part: SET2005150	20"W x 20"D ..... No Charge Part: SET2005154	

Contoured seating section continued on following pages.

**11. CONTOURED SEATING CONTINUED**

**11.2b. Quantum Captain Seat with Solid Seat Pan**

16"W x 16"D ..... No Charge Part: SET2005394	18"W x 20"D ..... No Charge Part: SET2005399	20"W x 22"D ..... No Charge Part: SET2005911
16"W x 18"D ..... No Charge Part: SET2005395	18"W x 22"D ..... No Charge Part: SET2005910	22"W x 18"D ..... No Charge Part: SET2005403
16"W x 20"D ..... No Charge Part: SET2005396	20"W x 16"D ..... No Charge Part: SET2005400	22"W x 20"D ..... No Charge Part: SET2005404
18"W x 16"D ..... No Charge Part: SET2005397	20"W x 18"D ..... No Charge Part: SET2005401	22"W x 22"D ..... No Charge Part: SET2005405
18"W x 18"D ..... No Charge Part: SET2005398	20"W x 20"D ..... No Charge Part: SET2005402	

**11.3. Seat to Floor Height**

Power adjustable seat height option is fixed at lowest seat to floor height.

**Lowest Seat to Floor**  
Part: PTOINDV3558

**Middle Seat to Floor**  
Part: PTOINDV3559

**Highest Seat to Floor**  
Part: PTOINDV3560

Low	Med	High
19.5"	20.5"	21.5"
Q Captain Seat to Floor Heights		

**11.4. Stealth Seat Cushion - Quick Selection**

**Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions can be found on the respective complete [Stealth Cushion Order Forms](#). This section allows a quick selection of a single cushion that matches the seat pan width and depth to alleviate listing all part numbers and codes on this power base form. See specific order form for non-matching sizes, alternate covers and liners, or additional covers and liners.**

**TRU-Comfort 2 SPP, Matching Seat Size**

Reference part: XRF000116. [See cushion order form](#) for part numbers, MSRP, HCPCS, and additional options. Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Part series: ST-TC2SPPXXYY. MSRP: \$441/\$551. HCPCS: E2607/E2608



**Essence SPP, Matching Seat Size**

Reference part: XRF000118. [See cushion order form](#) for part numbers, MSRP, HCPCS, and additional options. Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Part series: ST-ESSSPPXXYY. MSRP: \$441/\$551. HCPCS: E2607/E2608



**Select optional wedge kit below:**

Small Wedge Positioning Kit .....\$39  
Part: POP1909006-S  
14"W and smaller

Medium Wedge Positioning Kit.....\$72  
Part: POP1909006-M  
14"W and smaller

Large Wedge Positioning Kit..... \$105  
Part: POP1909006-L  
14"W and smaller

Small Wedge Positioning Kit .....\$39  
Part: POP1712340-S  
15"W and larger

Medium Wedge Positioning Kit.....\$72  
Part: POP1712340-M  
15"W and larger

Large Wedge Positioning Kit..... \$105  
Part: POP1712340-L  
15"W and larger

**Solution SPP, Matching Seat Size**

Reference part: XRF000113. [See cushion order form](#) for part numbers, MSRP, HCPCS, and additional options. Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner. Part series: ST-SOLSPPXXYY. MSRP: \$414/\$475. HCPCS: E2607/E2608



**Spectrum Foam SPP, Matching Seat Size**

Reference part: XRF000112. [See cushion order form](#) for part numbers, MSRP, HCPCS, and additional options. Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner. Part series: ST-SPCFSPXXYY. MSRP: \$507/\$585. HCPCS: E2607/E2608




**Spectrum Gel SPP, Matching Seat Size**

Reference part: XRF000117. [See cushion order form](#) for part numbers, MSRP, HCPCS, and additional options. Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner. Part series: ST-SPCGSPPXXYY. MSRP: \$507/\$585 HCPCS: E2607/E2608



**Premiere P, Matching Seat Size**

Reference part: XRF000122. [See cushion order form](#) for part numbers, MSRP, HCPCS, and additional options. Positioning with Coolcore, 4-way stretch, breathable cover. Part series: ST-PREMPXXYY. MSRP: \$358/\$524. HCPCS: E2605/E2606



**Glacial SP, Matching Seat Size**

Reference part: XRF000114. [See cushion order form](#) for part numbers, MSRP, HCPCS, and additional options. Skin protection with Coolcore, 4-way stretch, breathable cover. Part series: ST-GLSPXXYY. MSRP: \$364/\$485. HCPCS: E2603/E2604



**Zen SP, Matching Seat Size**

Reference part: XRF000115. [See cushion order form](#) for part numbers, MSRP, HCPCS, and additional options. Skin protection with Coolcore, 4-way stretch, breathable cover. Part series: ST-ZENSPXXYY. MSRP: \$348/\$468. HCPCS: E2603/E2604



**Simplicity G, Matching Seat Size**

Reference part: XRF000119. [See cushion order form](#) for part numbers, MSRP, HCPCS, and additional options. General use cushion with Stealth Tek 2, 2-way stretch, breathable cover. Part series: ST-SIMGXXYY. MSRP: \$122/\$182. HCPCS: E2601/E2602



Contoured seating section continued on following pages.

**11. CONTOURED SEATING CONTINUED**

**11.5. Armrests**

**REQUIRED.** Configured in pairs. Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.

**1a. Q-Captain, Flip-Up, Height Adjustable Armrests**

Includes full length, flat armpads.

Full Length, Left Armrest ..... \$90 ea.      Full Length, Right Armrest..... \$90 ea.  
 Part: MEC2201304. HCPCS: K0020      Part: MEC2201303. HCPCS: K0020

**1b. Drop-In, Height Adjustable, Armrests (T-Arm)**

Includes flat armpads.

Full Length, Drop-In Armrests..... \$79 ea.        Desk Length, Drop-In Armrests ..... \$79 ea.  
 L R Left part: XRF000298. Right part: XRF000299. HCPCS: E0973      L R Left part: XRF000300. Right part: XRF000301. HCPCS: E0973

**1c. Single Post, Height Adjustable, Removable Armrests**

Height adjustment from 11"-17" in half-inch increments.

Full Length, Single Post Armrests ..... \$79 ea.        Desk Length, Single Post Armrests ..... \$79 ea.  
 L R Left part: ARM2000906. Right part: ARM2000907. HCPCS: E0973      L R Left part: ARM2000908. Right part: ARM2000909. HCPCS: E0973

**1d. Bariatric, Height Adjustable, Removable Armrests**

Select armrests and armpads. Armrests are sold in same-size pairs.

Full Length, Bariatric Armrests, 8"-12" H ..... \$90 ea.        Full Length, Bariatric Armrests, 10"-14"H ..... \$90 ea.  
 L R Left part: FRMASMB6728. Right part: FRMASMB6729. HCPCS: E0973      L R Left part: FRMASMB6730. Right part: FRMASMB6731. HCPCS: E0973

**1e. Heavy Duty, 8"-12" Height Adjustable, Removable Armrests**

Select armrests and armpads. Armrests are sold in same-size pairs.

Full Length, Heavy Duty Armrests, 8"-12" H ..... \$90 ea.        Desk Length, Heavy Duty Armrests, 8"-12"H..... \$90 ea.  
 L R Left part: XRF000253. Right part: XRF000254. HCPCS: E0973      L R Left Part: XRF000255. Right part: XRF000256. HCPCS: E0973

**1f. Heavy Duty, 10"-14" Height Adjustable, Removable Armrests**

Select armrests and armpads. Armrests are sold in same-size pairs.

Full Length, Heavy Duty Armrests, 10"-14"H..... \$90 ea.        Desk Length, Heavy Duty Armrests, 10"-14" H ..... \$90 ea.  
 L R Left Part: XRF000257. Right part: XRF000258. HCPCS: E0973      L R Left part: XRF000259. Right part: XRF000260. HCPCS: E0973

**1g. Pediatric, Removable, Height Adjustable Armrests**

Height adjustment from 7.75" to 8.25". Includes desk length, straight armpads.

Pediatric, Removable Armrests ..... \$79 ea.  
 L R Left Part: ARM142876. Right part: ARM142877. HCPCS: E0973

**11.6. Armrest Height**

**REQUIRED with Bariatric, Heavy Duty, and Pediatric armrests.**

**2. Left Armrest Height**

7.75" H PTOINDV3722	8.25" H PTOINDV3724	10" H PTOINDV3537	12" H PTOINDV3539	14" H PTOINDV3541	16" H PTOINDV3543
8" H PTOINDV3535	9" H PTOINDV3536	11" H PTOINDV3538	13" H PTOINDV3540	15" H PTOINDV3542	

**3. Right Armrest Height**

7.75" H PTOINDV3721	8.25" H PTOINDV3723	10" H PTOINDV3452	12" H PTOINDV3454	14" H PTOINDV3456	16" H PTOINDV3458
8" H PTOINDV3450	9" H PTOINDV3451	11" H PTOINDV3453	13" H PTOINDV3455	15" H PTOINDV3457	

**11.7. Armpads**

**REQUIRED except with Flip-Up, Drop-In, and Pediatric armrests.**

**1a. Straight Armpads for Single Post Armrests**

Straight, Full Length Left and Right ..... No Charge      Straight, Desk Length Left, Full Length Right ..... No Charge  
 Part: ARM134275      Part: ARM134280

Straight, Full Length Left, Desk Length Right ..... No Charge      Straight, Desk Length Left and Right..... No Charge  
 Part: ARM134281      Part: ARM134254

**1b. Waterfall Armpads for Single Post Armrests**

Waterfall, Full Length Left and Right..... No Charge      Waterfall, Desk Length Left, Full Length Right ..... No Charge  
 Part: ARM134285      Part: ARM134288

Waterfall, Full Length Left, Desk Length Right ..... No Charge      Waterfall, Desk Length Left and Right ..... No Charge  
 Part: ARM134289      Part: ARM134282

**1c. Falcon Gel Armpads for Single Post Armrests**

Falcon Gel, Full Length Left and Right ..... \$311.22      Falcon Gel, Desk Length Left, Full Length Right..... \$287.93  
 Part: ARM136438      Part: ARM136439

Falcon Gel, Full Length Left, Desk Length Right..... \$287.93      Falcon Gel, Desk Length Left and Right ..... \$265.02  
 Part: ARM136442      Part: ARM136437

Contoured seating section continued on following pages.

**11. CONTOURED SEATING CONTINUED**

**1d. Armpads for Bariatric or Heavy Duty Armrests**

Bariatric Armpads ..... No Charge  
 Part: UPSASMB1027  
 Only for use with Bariatric Armrests.

Straight Armpads..... No Charge  
 Part: UPSASMB1021  
 Only for use with Heavy Duty Armrests

Waterfall Armpads ..... No Charge  
 Part: UPSASMB1024  
 Only for use with Heavy Duty Armrests

**11.8. Quad Release Lever**

Optional. Available with Single Post and Pediatric armrests.

Quad Release Lever, Pair.....No Charge  
 Part: FRMASMB10517

**12. LAP BELTS**

**REQUIRED.** Select a lap belt below.

**12.1a. Standard Lap Belts**

38" Lap Belt Part: ACCBELT1006	50" Lap Belt Part: ACCBELT1000	60" Lap Belt Part: ACCBELT1002	70" Lap Belt Part: ACCBELT1003	80" Lap Belt Part: ACCBELT1004	90" Lap Belt Part: ACCBELT1001
-----------------------------------	-----------------------------------	-----------------------------------	-----------------------------------	-----------------------------------	-----------------------------------

**12.1b. Optional Lap Belts**

Stealth Positioning 60" Single Pull, Push Button, Medium..... \$50 (7.5") Padded, 1.5" Belt (IFPB5213) Part: ST-IFPB5213-Q	Stealth Positioning 60" Single Pull, Push Button, Extra Large..... \$50 Padded (11.5"), 2" Belt (IFPB5315) Part: ST-IFPB5315-Q
Stealth Positioning 60" Single Pull, Push Button, Large ..... \$50 Padded (9.5"), 2" Belt (IFPB5314) Part: ST-IFPB5314-Q	60" Sure-Lok Lap Belt (Non-retractable)..... \$120 Required with ISO 7176-19 Compliant Transit Option. Coming soon. Part: ACCBELT1016

**13. LEGRESTS**

Optional. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.

**13.1a. Basic Center Mount Foot Platform**

Select legrest length in following section.

**1. Center Mount Foot Platform Type**

Basic Center Mount.....No Charge  
 Part: XRF000323

**2. Basic Center Mount Foot Platform Extension and Angle Adjustable Footplates**

See legrest length matrix for compatibility.

Short Extension, 5"x9" Footplates .....No Charge Part: MEC1906675	Long Extension, 5"x9" Footplates .....No Charge Part: MEC1906677
Medium Extension, 5"x9" Footplates .....No Charge Part: MEC1906680	

**3. Basic Center Mount Calf Pads**

Calf pads are standard with a power articulating foot platform.

Small Calf Pads, Basic Center Mount..... \$158 Part: RIG2104369	Large Calf Pads, Basic Center Mount ..... \$158 Part: RIG2202997 Only available with Long Length Extension.
--	---

**4. Basic Center Mount Additional Accessories**

Basic Center Mount Footplate Connector ..... \$16  
 Part: RIG2104240

**13.1b. Heavy Duty, Drop-In Legrests**

Select a configuration, a legrest, extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

**1. Drop-In Legrest Selection**

<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 70° HD Drop-in for 8.5"-12" Lower Extension ..... No Charge Pair part: FRMASMB4762. Left: RIG137706. Right: RIG137707	<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 70° HD Drop-in 13"-18" or 16"-21" Lower ..... No Charge Extension Pair part: FRMASMB4853. Left: RIG137708. Right: RIG137709
---	--

**2. Lower Extension**

<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length Pair part: XRF000101. Left: XRF000102. Right: XRF000103.	<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 16"-21" HD Drop-Ins/18.5"-22.5" ELR, Ext Length Pair part: FRMASMB4754. L: FRMASMB4725. R: FRMASMB5393
<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R 13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length Pair part: XRF000104. Left: XRF000105. Right: XRF000106	

Legrests section continued on following page.

**13. LEGRESTS CONTINUED**

**3. Angle Adjustable Footplate(s) with Heel Loops**

- P L R** 6"Wx8"D<sup>1</sup> ..... \$90 ea.  
*Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795.*  
*HCPCS: K0040+E0951*
- P L R** 6"Wx11"D<sup>1</sup> ..... \$90 ea.  
*Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797.*  
*HCPCS: K0040+E0951*  
*May cause interference on seat widths less than 18".*

- P L R** 8"Wx11"D<sup>2</sup> ..... \$90 ea.  
*Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799.*  
*HCPCS: K0040+E0951*  
*May cause interference on seat widths less than 20".*

**13.2. Legrest Accessories**

- Heel Loops (pair) ..... \$53  
*Part: FRMASMB7873. HCPCS: E0951*  
*Standard on angle adjustable and multi axis footplates when noted.*
- Wheel Bumper ..... \$37  
*Part: XRF000100*

**13.3. Set Legrest Extension Length**

**REQUIRED with matching left and right legrest configuration. Part number is PTOINDV plus number following # below.** All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference. Non-matching or individual legrest configurations lengths are not set by the manufacturer.

**1a. Centermount Extension Length**

2" L (#3989)	5" L (#3992)	8" L (#3680)	9.5" L (#3683)	11" L (#3686)	12.5" L (#3689)	14" L (#3692)	15.5" L (#3695)	17" L (#3698)	18.5" L (#3701)
3" L (#3990)	6" L (#3993)	8.5" L (#3681)	10" L (#3684)	11.5" L (#3687)	13" L (#3690)	14.5" L (#3693)	16" L (#3696)	17.5" L (#3699)	19" L (#3702)
4" L (#3991)	7" L (#3994)	9" L (#3682)	10.5" L (#3685)	12" L (#3688)	13.5" L (#3691)	15" L (#3694)	16.5" L (#3697)	18" L (#3700)	

**1b. Front Rigging Extension Length**

Non-matching or individual legrest configurations lengths are not set by the manufacturer. See previous section for center mount extension lengths.

8" L (#3472)	9" L (#3474)	10" L (#3476)	11" L (#3478)	12" L (#3480)	13" L (#3482)	14" L (#3484)	15" L (#3486)	16" L (#3488)	17" L (#3490)
8.5" L (#3473)	9.5" L (#3475)	10.5" L (#3477)	11.5" L (#3479)	12.5" L (#3481)	13.5" L (#3483)	14.5" L (#3485)	15.5" L (#3487)	16.5" L (#3489)	

**14. POWER LIFT ELECTRONICS**

**REQUIRED with 12" Power Lift for Captain Seat with iLevel.**

- Single Actuator Through Controller and Interface ..... \$2,163  
*Part: XRF000335. HCPCS: E2310*

**15. ACCESSORIES**

*Optional.*

- |   |  |
|---|--|
| <ul style="list-style-type: none"> <li>Cup Holder, Q-Captain ..... No Charge<br/> <i>Part: ACC2006355</i></li> <li>Front Mounted Cell Phone/ Media Holder, X-Grip, Q-Captain..... \$74<br/> <i>Part: ACC2103691</i></li> <li>Rear View Mirror, Q-Captain ..... \$69<br/> <i>Part: ACC2103692</i></li> <li>Oxygen Holder.....\$118<br/> <i>Part: ACC2006353. HCPCS: E2208</i></li> <li>Cane/Crutch Holder .....\$110<br/> <i>Part: ACC2006352. HCPCS: E2207</i></li> <li>Rear Basket .....\$118<br/> <i>Part: ACC2006351</i></li> <li>Saddle Bag.....\$40<br/> <i>Part: DWR1020D010</i></li> <li>Large, Black, Weather Cover ..... \$202<br/> <i>Part: ACCCOVR1021</i></li> <li>Quantum Backpack ..... \$168<br/> <i>Part: ACC1905853</i></li> <li>Grip Lap Board (20005)..... \$38.95<br/> <i>Part: ACC152992</i></li> <li>Grip Drink Holder (20085) ..... \$15.50<br/> <i>Part: ACC152993</i></li> </ul> | <ul style="list-style-type: none"> <li>Reflector Kit ..... \$53<br/> <i>Part: ACC124670</i></li> <li>Clothing Guard (Pair)..... \$237<br/> <i>Part: ACC125007</i><br/> <i>Available with Flip-Up armrests only. Not available with Elevating Legrests.</i></li> <li>Transfer Handles .....\$105 ea.<br/> <i>Part: MEC1809930. HCPCS: E0705</i><br/> <i>Available with Flip-Up armrests only. Not available with Elevating Legrests.</i></li> <li>Privacy Bag.....\$84<br/> <i>Part: ACC2005723</i></li> <li>Personal Item Pouch, Rail Mounted.....\$84<br/> <i>Part: ACC2101017</i></li> <li>Personal Item Pouch, Arm Mounted .....\$84<br/> <i>Part: ACC2103694</i></li> <li>Glove Box, Rail Mounted ..... \$158<br/> <i>Part: ACC125011</i></li> <li>Glove Box, Arm Mounted..... \$158<br/> <i>Part: ACC2103693</i></li> <li>Q-Logic EX/NE Programmer ..... \$573<br/> <i>Part: ELEASMB6976</i></li> <li>Q-Logic PC Programming Disk with Cable ..... \$445<br/> <i>Part: ELEASMB5216</i></li> </ul> |
|---|--|

**16. BRANCH PATCHES**

Optional.

Air Force Patch ..... \$21 Part: ACC140240	Marines Patch..... \$21 Part: ACC140243
Army Patch ..... \$21 Part: ACC140241	Navy Patch ..... \$21 Part: ACC140244
Coast Guard Patch ..... \$21 Part: ACC140239	Purple Heart Patch ..... \$21 Part: ACC140238



**17. POSITIONING COMPONENTS**

Optional. More options are available on the [Positioning Components Order Form](#). Please provide quantity next to check box where applicable.

**17.1. Stealth Unilink Pelvic/Thigh Guide**

**1. UniLink Pelvic/Thigh Hardware**

UniLink, 5" Height Adjustable, ..... \$265 ea. Removable Mounting Hardware Part: ST-UTTP5-Q. HCPCS: E1028	Flip Down with LEFT Full Surface Contact ..... \$281 ea. Feature (TWBADD-L & LPHW407) Part: POS140428. HCPCS: E1028
UniLink, 7" Height Adjustable, ..... \$276 ea. Removable Mounting Hardware Part: ST-UTTP7-Q. HCPCS: E1028	Flip Down with RIGHT Full Surface Contact..... \$281 ea. Feature (TWBADD-R & LPHW407) Part: POS140429. HCPCS: E1028

**2. UniLink Pelvic/Thigh Pads**

^Please verify the appropriate billing code to be used for this item with all third-party payors.

UniLink, 4" x 4" Pad ..... \$106 ea. Part: ST-UTP544056100-Q. HCPCS: E0953/E0956^	UniLink, 4" x 12" Pad ..... \$106 ea. Part: ST-UTP54D056100-Q. HCPCS: E0953/E0956^
UniLink, 4" x 8" Pad..... \$106 ea. Part: ST-UTP548056100-Q. HCPCS: E0953/E0956^	

**3. UniLink Pelvic/Thigh Additional Links**

UniLink, 1.5" Link ..... \$81 ea. Part: ST-UIL001150A000-Q	UniLink, 2" Link..... \$81 ea. Part: ST-UIL001200A000-Q
---	--

**4. UniLink Pelvic/Thigh Guide Installation**

Install on Left Part: PTOINDV3585	Install on Right Part: PTOINDV3586	Do Not Install Part: PTOINDV3719
--------------------------------------	---------------------------------------	-------------------------------------

**17.2. Stealth Pelvic/Thigh Guides**

Pelvic/thigh guides may interfere with the joystick receiver block and down posts.

**1. Pelvic/Thigh Guide Bracket**

Flip Down Bracket (TWBTLADD-HW) ..... \$281 ea. Part: INDPART2804. HCPCS: E1028	Removable Bracket with 1" Offset ..... \$382 ea. (LP430 & LHW-121) Part: INDPART2854. HCPCS: E1028
Removable Bracket (LP430)..... \$284 ea. Part: INDPART2805. HCPCS: E1028	

**2. Pelvic/Thigh Guide Adapter**

2" Height Adapter Plate (LHW-122)..... \$28 Part: INDPART2821. HCPCS: K0108 Height adapter is needed with TRU-Comfort seat cushions. Not available with the flip down bracket.	2 1/4" Adapter Plate (LHW-129) ..... \$57 Part: POS140430. HCPCS: K0108 Only available with the flip down full surface bracket.
--	---

**3. Pelvic/Thigh Guide Pad**

^Please verify the appropriate billing code to be used for this item with all third-party payors.

Short Pad, 4"x4" (TWBADD-4X4)..... \$81 ea. Part: INDPART2801. HCPCS: E0953/E0956^	Long Pad, 4"x12" (TWBADD-4X12)..... \$81 ea. Part: INDPART2803. HCPCS: E0953/E0956^
Medium Pad, 4"x8" (TWBADD-4X8)..... \$81 ea. Part: INDPART2802. HCPCS: E0953/E0956^	

**4. Pelvic/Thigh Guide Installation**

Install on Left Part: PTOINDV3585	Install on Right Part: PTOINDV3586	Do Not Install Part: PTOINDV3719
--------------------------------------	---------------------------------------	-------------------------------------

**17.3. Stealth Elbow Support**

**1. Elbow Support Bracket**

Left, Swing-Away, Elbow Support Bracket, Tube ..... \$199 ea. Style Arms (ARMS-ES-L) Part: POS124064. HCPCS: E1028	Right, Swing-Away, Elbow Support Bracket, Tube ..... \$199 ea. Style Arms (ARMS-ES-R) Part: POS124066. HCPCS: E1028
--	---

**2. Elbow Support Pad**

3"Wx3"H Elbow Pad (ARMS-ES-3X3) ..... \$57 ea. Part: POS124065. HCPCS: K0108	4"Wx4"H Elbow Pad (ARMS-ES-4X4)..... \$57 ea. Part: POS124067. HCPCS: K0108
---	--

Positioning components section continued on following pages.



**17. POSITIONING COMPONENTS CONTINUED**

**17.4. Stealth Foot Boxes**

Stealth Small Foot Boxes, 5"Wx8"Dx7"H..... \$539  
 (SPLE107) (pair)  
 Part: INDPART2786. HCPCS: E0954

Stealth Large Foot Boxes, 7"Wx12"Dx11"H..... \$539  
 (SPLE109) (pair)  
 Part: INDPART2788. HCPCS: E0954

Stealth Medium Foot Boxes, 6"Wx10"Dx9"H ..... \$539  
 (SPLE108) (pair)  
 Part: INDPART2787. HCPCS: E0954

**17.5. Stealth Positioning Shoe Holder**

Small Shoe Holder, 3/4" Webbing (IFSH4211)..... \$87 ea.  
 Part: ST-IFSH4211-Q. HCPCS: E0951+E0952

Large Shoe Holder, 3/4" Webbing (IFSH4213) ..... \$87 ea.  
 Part: ST-IFSH4213-Q. HCPCS: E0951+E0952

Medium Shoe Holder, 3/4" Webbing (IFSH4212) ..... \$87 ea.  
 Part: ST-IFSH4212-Q. HCPCS: E0951+E0952

**18. LAP TRAY**

**18.1. Simple Solutions Lap Tray**

Mounting bracket is required with lap tray selection. See chart below for compatibility.

**1. Lap Tray**

Simple Solutions Small Tray with Cut Out ..... \$309  
 Part: ST-SSLXS010001171-Q. HCPCS: E0950

Simple Solutions Large Tray with Cut Out..... \$309  
 Part: ST-SSLXL010001171-Q. HCPCS: E0950

Simple Solutions Medium Tray with Cut Out..... \$309  
 Part: ST-SSLXM010001171-Q. HCPCS: E0950

**2. Mounting Hardware**

Swing Away Mounting Hardware..... \$302  
 Part: ST-SSLXHD10000060-Q

Simple Solutions Lap Tray Armrest Compatibility	Width											
	10	11	12	13	14	16	17	18	19	20	21	22
Flip-Back	S	S	S	S	S	S	S	M	ML	ML	L	L
Flip-Back, Cane Mounted	S	S	S	S	S	S	M	ML	ML	L	L	
Drop-In, T-Arm	S	S	S	S	S	S	M	ML	ML	L	L	
Two-Post, Flip-Back	S	S	S	S	S	S	M	ML	ML	L	L	
Cantilever	S	S	S	S	S	M	M	ML	ML	L	L	
Q-Captain, Flip-Up						S	M	ML		L		
Single Post/ Heavy Duty/ Bariatric	S	S	S	S	S	S	M	S	M	L	L	L

S:ST-SSLXS010001171-Q    M:ST-SSLXM010001171-Q    L:ST-SSLXL010001171-Q  
 :Compatible    :Armrest configuration may not hold tray position  
 \*Some armrests listed may not be available on this order form.