



Quantum Outback Order Form with TRU Balance® 3 Positioning

Account Number: _____ Date: _____
 Provider Name: _____
 Contact: _____
 Phone: _____ Fax: _____
 Email: _____
 PO Number: _____
 Quote Name/Marked for: _____
 Ship to Address: _____
 City: _____ State: _____ Zip: _____



INTRODUCTION

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.


This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email (quantumorders@pridemobility.com). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering.

Thank you for choosing Quantum!

PATIENT INFORMATION

Weight: _____ Height: _____ A. Shoulder Width: _____ D. Max Sitting Height: _____ G. Thigh Depth: _____
 B. Chest Width: _____ E. Shoulder Height: _____ H. Lower Leg Length: _____
 C. Hip Width: _____ F. Axilla Height: _____ I. Elbow Height: _____



1. BASE MODEL

REQUIRED. 300 lb. weight capacity with TRU-Balance 3 positioning, 14" pneumatic drive wheels, and 8 amp off board charger. Dimensions of base without legrests: 28.82"Wx42"L.

Quantum Outback \$19,595
 Part: Q OUTBACK. HCPCS: K0899 - Power mobility device, not coded by PDAC

2. BASE COLOR

REQUIRED. Select one shroud color.

Cadet Base Shroud Color No Charge Part: DGN2301319-TGY3-PCL2	Desert Camouflage Base Shroud Color <i>Coming soon</i> Part: DGN2301319-WCM11-PCL2
Tan Base Shroud Color No Charge Part: DGN2301319-TTN1-PCL2	Tropic Lightning Camouflage Base Shroud Color <i>Coming soon</i> Part: DGN2301319-WCM12-PCL2
Black Camouflage Base Shroud Color <i>Coming soon</i> Part: DGN2301319-WCM9-PCL2	Jungle Camouflage Base Shroud Color <i>Coming soon</i> Part: DGN2301319-WCM13-PCL2
White Camouflage Base Shroud Color <i>Coming soon</i> Part: DGN2301319-WCM10-PCL2	

3. MOTORS

REQUIRED.

5 MPH Motors No Charge Part: MOT2205971	6 MPH Motors No Charge <i>Coming soon.</i> Part: MOT2300107
--	--

4. BASE MOUNTED LIGHTS

Optional.

Omit Base Mounted Lights No Charge Part: DGN2400803
--

5. QUANTUM BACKUP CAMERA

Optional. The backup camera may interfere with armrest mounted accessories such as media device holder.

Quantum Backup Camera \$368 Part: MEC1910365

1. Optional Camera Mounting Bracket

Backup Camera Headrest Bracket \$21 Part: MEC1910373 Works with .5" x .5" vertical supports only.	Backup Camera Track Mount Bracket \$27 Part: MEC2104227
---	--

2. Optional Display Mounting Bracket

Backup Camera Q-Logic Joystick Mount Bracket \$105 Part: MEC2106123
--



6. BATTERIES

Optional. The powerchair requires 2 batteries and will be installed if ordered.

Group 27, Interceptor 2, Gel Battery \$519 ea. Part: BAT2401900
--

7. ELECTRONICS

REQUIRED. Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Form](#).

7.1a. Omit Joystick

No Joystick, Setup for Q-Logic Electronics Part: XRF000010 Requires selection of Enhanced Display

7.1b. Q-Logic 3e Joystick

Q-Logic 3e Joystick No Charge Part: CTL1910600 Can operate full light package. Not available with alternative or specialty controls.
--

7.1c. Q-Logic 3 EX Joystick

Q-Logic 3 expandable capacity controls allow use of specialty controls and multiple actuators.

1. Q-Logic 3 Required Components

Q-Logic 3 controls require the controller module and expandable harness below. If a Q-Logic joystick or specialty control is selected, Customer Service will quote the order with these items.

Q-Logic 3 EX Controller Module \$730 Part: XRF000002. HCPCS: E2377	Q-Logic 3 Expandable Controller Harness \$458 Part: XRF000264. HCPCS: E2313
---	--

2. Q-Logic 3 Joystick

Q-Logic 3 EX Joystick \$840 Part: CTL167724	Q-Logic 3 EX Joystick for Full Light Package \$840 Part: CTL167725 Required with full light packages.
--	---



8. JOYSTICK HANDLES

Optional. The standard Q-Logic 3 joystick knob can be seen in an image in the previous section. An optional handle selected below will be mounted to the joystick and the standard handle will be shipped with the unit.

- | | |
|--|--|
| Mushroom Handle (B) \$90
Part: ACCASMB1860. HCPCS: E2323 | Extended I-Handle (G) \$90
Part: ACCASMB1857. HCPCS: E2323 |
| Chin Cup (C) \$90
Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control) | Bodypoint 3" Goal Post Handle \$104
Part: PLSKNCP1057. HCPCS: E2323 |
| T-Bar (D) \$126
Part: ACCASMB1858. HCPCS: E2323 | Bodypoint 4" Goal Post Handle \$104
Part: PLSKNCP1058. HCPCS: E2323 |
| Medium Ball Handle (E) \$69
Part: ACCASMB1859. HCPCS: E2323 | Q-Logic 3 Carrot Handle (A) No Charge
Part: DGN124308 |
| Soft Ball Handle (F) \$80
Part: ACCASMB1856. HCPCS: E2323 | |



9. JOYSTICK MOUNTING BRACKETS

Optional. Select a mount to be used with a joystick selected in previous section.

9.1.a. Fixed, In-line Joystick Bracket

- | |
|--|
| Fixed, In-line, Bracket Standard
Part: FRMASMB10617 |
|--|




9.1.b. Swing-Away (Retractable) Joystick Bracket

Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation.

- | |
|---|
| Swing-away, Q-Logic 3 \$319
Part: ACC1908434. HCPCS: E1028 |
|---|

9.1.c. Basic Swing Away Joystick Mount

- | | |
|--|---|
| Basic Swing Away, Q-Logic, Left Side \$245
Part: ACC2201041. HCPCS: E1028 | Basic Swing Away, Q-Logic, Right Side \$245
Part: ACC2201042. HCPCS: E1028 |
|--|---|



9.1.d. Specialty Joystick Mounts

Select a bracket and side to mount it.

- | | |
|--|--|
| Retract4 Mount \$368
Part: XRF000003. HCPCS: E1028
Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted. | Stealth Swing-Away/Flip Down Bracket \$263
Part: XRF000005. HCPCS: E1028
Allows the joystick to pivot to side of armrest, or flip-down if bracket is adjusted. |
| Multi-Axis Mount \$126
Part: XRF000004
Allows height adjustment of joystick. | |



9.2. Joystick Mount Location

- | | |
|--------------------------------------|---------------------------------------|
| Left Side Mount
Part: PTOINDV1225 | Right Side Mount
Part: PTOINDV1226 |
|--------------------------------------|---------------------------------------|

10. TRU BALANCE 3 SEATING

10.1. TRU Balance 3 Configuration

REQUIRED. Select one seating system. Be sure to complete the seat dimensions, seat to floor, and back selections.

- | | |
|---|---|
| TRU Balance 3 Power Tilt Only \$8,202
Part: PTOINDV3443. HCPCS: E1002
50° Center of Gravity (CG) power tilt seating system. | TRU Balance 3 Power Tilt & Recline \$12,615
Part: PTOINDV3446. HCPCS: E1007
50° CG power tilt with 175° of power recline. Full drive lockout at 20°. |
| TRU Balance 3 Power Recline Only \$8,234
Part: PTOINDV3444. HCPCS: E1004
85° to 175° of power recline. | TRU Balance 3 Static Seat \$2,410
Part: PTOINDV3442
Must select power lift. 85° to 105° of limited recline (-5/ +15) via gas strut within back cover. |

10.2. 10" Power Lift

Optional. Select one option below. ^Please verify the appropriate billing code to be used for this item with all third-party payors

- | |
|--|
| TRU Balance 3 Power HD Lift \$3,725
Part: PTOINDV3797. HCPCS: E2298/E2300^
10" of scissor lift. 22"W x 22"D maximum seat size. |
|--|

10.3. Occupied Transit

Optional. Available in seat dimensions of 16" to 20" wide and 16"-20" deep. Requires TRU Balance 3 armrests and 60" Sure-lok lap belt. Only available with swing-away legrests and center mount foot platforms.

- | | |
|--|---|
| Occupied Transit \$289
Coming soon. Part: KIT152047 | To determine if your chair is WC-19 compliant (occupied), please refer to this website: https://wc-transportation-safety.umtri.umich.edu/ |
|--|---|

10.4. Seat Width

REQUIRED with seat selection. Minimum seat width is 16". Maximum seat width is 22".

- | | | | | | | |
|------------------|------------------|------------------|------------------|------------------|------------------|------------------|
| 16" W
(#1006) | 17" W
(#1007) | 18" W
(#1008) | 19" W
(#1009) | 20" W
(#1010) | 21" W
(#1011) | 22" W
(#1012) |
|------------------|------------------|------------------|------------------|------------------|------------------|------------------|

TRU Balance 3 seating section continued on following pages.

10. TRU BALANCE 3 SEATING CONTINUED

10.5. Seat Width Package

REQUIRED with seat selection.

16"-20" Width Package
Part: XRF000008

20"-24" Width Package
Part: XRF000158

10.6. Seat Depth

REQUIRED with seat selection. Minimum seat depth is 16". Maximum seat depth is 22".

16" D (#1006)	17" D (#1007)	18" D (#1008)	19" D (#1009)	20" D (#1010)	21" D (#1011)	22" D (#1012)
------------------	------------------	------------------	------------------	------------------	------------------	------------------

10.7. Seat to Floor Height

REQUIRED.

Lowest Seat to Floor
Part: PTOINDV3558

Middle Seat to Floor
Part: PTOINDV3559

Highest Seat to Floor
Part: PTOINDV3560

Low	Med	High
19.5"	20.5"	21.5"
Q Captain Seat to Floor Heights		

10.8. Stealth Seat Cushion - Quick Selection

Optional. Specific cushion/component sizes available, part numbers, weight capacities, MSRP, HCPCS code, and special order cushions can be found on the respective complete [Stealth Cushion Order Forms](#). This section allows a quick selection of a single cushion that matches the seat pan width and depth to alleviate listing all part numbers and codes on this power base form. See specific order form for non-matching sizes, alternate covers and liners, or additional covers and liners.

TRU-Comfort 2 SPP, Matching Seat Size

Reference part: XRF000116. [See cushion order form](#) for part numbers, MSRP, HCPCS, and additional options.
Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.
Part series: ST-TC2SPPXXYY. MSRP: \$441/\$551. HCPCS: E2607/E2608



Essence SPP, Matching Seat Size

Reference part: XRF000118. [See cushion order form](#) for part numbers, MSRP, HCPCS, and additional options.
Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.
Part series: ST-ESSSPPXXYY. MSRP: \$441/\$551. HCPCS: E2607/E2608



Select optional wedge kit below:

Small Wedge Positioning Kit\$39
Part: POP1909006-S
14"W and smaller

Medium Wedge Positioning Kit.....\$72
Part: POP1909006-M
14"W and smaller

Large Wedge Positioning Kit..... \$105
Part: POP1909006-L
14"W and smaller

Small Wedge Positioning Kit\$39
Part: POP1712340-S
15"W and larger

Medium Wedge Positioning Kit.....\$72
Part: POP1712340-M
15"W and larger

Large Wedge Positioning Kit..... \$105
Part: POP1712340-L
15"W and larger

Solution SPP, Matching Seat Size

Reference part: XRF000113. [See cushion order form](#) for part numbers, MSRP, HCPCS, and additional options.
Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.
Part series: ST-SOLSPPXXYY. MSRP: \$414/\$475. HCPCS: E2607/E2608



Spectrum Foam SPP, Matching Seat Size

Reference part: XRF000112. [See cushion order form](#) for part numbers, MSRP, HCPCS, and additional options.
Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.
Part series: ST-SPCFSPXXYY. MSRP: \$507/\$585. HCPCS: E2607/E2608




Spectrum Gel SPP, Matching Seat Size

Reference part: XRF000117. [See cushion order form](#) for part numbers, MSRP, HCPCS, and additional options.
Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.
Part series: ST-SPCGSPPXXYY. MSRP: \$507/\$585 HCPCS: E2607/E2608



Premiere P, Matching Seat Size

Reference part: XRF000122. [See cushion order form](#) for part numbers, MSRP, HCPCS, and additional options.
Positioning with Coolcore, 4-way stretch, breathable cover.
Part series: ST-PREMPXXYY. MSRP: \$358/\$524. HCPCS: E2605/E2606




Glacial SP, Matching Seat Size

Reference part: XRF000114. [See cushion order form](#) for part numbers, MSRP, HCPCS, and additional options.
Skin protection with Coolcore, 4-way stretch, breathable cover.
Part series: ST-GLSPXXYY. MSRP: \$364/\$485. HCPCS: E2603/E2604



Zen SP, Matching Seat Size

Reference part: XRF000115. [See cushion order form](#) for part numbers, MSRP, HCPCS, and additional options.
Skin protection with Coolcore, 4-way stretch, breathable cover.
Part series: ST-ZENSPPXXYY. MSRP: \$348/\$468. HCPCS: E2603/E2604



Simplicity G, Matching Seat Size

Reference part: XRF000111. [See cushion order form](#) for part numbers, MSRP, HCPCS, and additional options.
General use cushion with Stealth Tek 2, 2-way stretch, breathable cover.
Part series: ST-SIMGXXYY. MSRP: \$122/\$182. HCPCS: E2601/E2602



TRU Balance 3 seating section continued on following pages.

10. TRU BALANCE 3 SEATING CONTINUED

10.9. Back Type Selection

REQUIRED. Match a choice in this section with the options in the following back cushion configuration section.

Redesigned, Static Type Back

Part: XRF000014

12" and 24" back width will use the original TRU Balance 3 Back assembly. Select a back cushion configuration in the subsequent section.

Back Cane Option

Part: XRF000015

Available only with TRU Balance 3 Static Seat and Power Tilt (1 actuator) selections. Select cane type, back cane height, and back cushion configuration below.

Redesigned, Flip Forward Back \$184

Part: XRF000016

12" and 24" back width will use the original TRU Balance 3 Back assembly. Not available with recline systems. Only available with 2-post, flip back or pediatric armrests. Not available with 12"-14" width. Select a back cushion type in the subsequent section.

10.10. Back Cane Selection

REQUIRED when Back Cane Option is selected above. Select type and height of back canes.

1. Back Cane Type

Angle Adjustable Back Cane No Charge

Part: MEC140266

2" Endomorph (1" per side) Back Cane \$74

Part: XRF000269

Flip Forward Back Cane \$53

Part: XRF000273

4" Endomorph (2" per side) Back Cane \$74

Part: XRF000271

2. Back Cane Height

16"H Cane

Part: FRMASMB15834

20"H Cane

Part: FRMASMB15838

24"H Cane

Part: FRMASMB15842

17"H Cane

Part: FRMASMB15835

21"H Cane

Part: FRMASMB15839

25"H Cane

Part: FRMASMB15843

10.11. Back Cushion Configuration

REQUIRED. Sport Back and TRU-Comfort Back selections up to 22" back width will be the Sport 4 or TRU-Comfort 4 back cushions. Sport back and TRU-Comfort back selections of 24" back width will be the Sport or TRU-Comfort 2 back cushions. Select a configuration type and applicable dimensions. See the Sport, Q-Back, and TRU-Comfort back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time.

1a. Sport Back/Omit/After-Market Interface

Sport Back No Charge

Part: XRF000017

Complete following width package and size.

After-market Back Interface \$535

Part: ST-SSAPT4-Q. Select a set width in following section, range of 12"-21". Overall seat depth may need to be adjusted to accommodate certain backs. May interfere with gel/waterfall arm pads and arm through. Does not include brackets to mount after-market backs.

Omit Back Pan/Cushion No Charge

Part: PTOINDV1287

Select a following width package.

1b. TRU-Comfort Back

TRU-Comfort Back (Static/Flip Forward), 12" Width \$744

Part: SET123622. HCPCS: E2620

Complete following width package and size.

TRU-Comfort Back (Static/Flip Forward), 22" Width \$855

Part: TC4V22W32W. HCPCS: E2621

Complete following width package and size.

TRU-Comfort Back (Static/Flip Forward), 14" to 20" Width \$744

Part: TC4V10W21W. HCPCS: E2620

Complete following width package and size.

TRU-Comfort Back, Back Cane Option, 22" Width \$855

Part: TC4V22W32WBC. HCPCS: E2621

Complete following width package and size.

TRU-Comfort Back, Back Cane Option, 16" to 20" Width \$744

Part: TC4V10W21WBC. HCPCS: E2620

Complete following width package and size.

1c. ADI Back, AL Series, 4-Point, Active Contour

Available with Back Cane Option only. Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Active Contour, 14" Wide x 13" or 14" Overall Height \$622

Part: ALBR14LA4

Active Contour, 18" Wide x 21" or 22" Overall Height \$622

Part: ALBR18TA4. HCPCS: E2613

Active Contour, 14" Wide x 16" or 17" Overall Height \$622

Part: ALBR14MA4

Active Contour, 18" Wide x 23" or 24" Overall Height \$622

Part: ALBR18XA4. HCPCS: E2613

Active Contour, 16" Wide x 19" or 20" Overall Height \$622

Part: ALBR16SA4. HCPCS: E2613

Active Contour, 20" Wide x 19" or 20" Overall Height \$622

Part: ALBR20SA4. HCPCS: E2613

Active Contour, 16" Wide x 21" or 22" Overall Height \$622

Part: ALBR16TA4. HCPCS: E2613

Active Contour, 20" Wide x 21" or 22" Overall Height \$622

Part: ALBR20TA4. HCPCS: E2613

Active Contour, 16" Wide x 23" or 24" Overall Height \$622

Part: ALBR16XA4. HCPCS: E2613

Active Contour, 20" Wide x 23" or 24" Overall Height \$622

Part: ALBR20XA4. HCPCS: E2613

Active Contour, 18" Wide x 19" or 20" Overall Height \$622

Part: ALBR18SA4. HCPCS: E2613

TRU Balance 3 seating section continued on following pages.

10. TRU BALANCE 3 SEATING CONTINUED

1d. ADI Back, AL Series, 4-Point, Deep Contour

Available with Back Cane Option only. Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Deep Contour, 14" Wide x 13" or 14" Overall Height..... \$735 Part: ALBR14LD4	Deep Contour, 18" Wide x 19" or 20" Overall Height..... \$735 Part: ALBR18SD4. HCPCS: E2615
Deep Contour, 14" Wide x 16" or 17" Overall Height..... \$735 Part: ALBR14MD4	Deep Contour, 18" Wide x 21" or 22" Overall Height..... \$735 Part: ALBR18TD4. HCPCS: E2615
Deep Contour, 16" Wide x 19" or 20" Overall Height..... \$735 Part: ALBR16SD4. HCPCS: E2615	Deep Contour, 20" Wide x 19" or 20" Overall Height \$735 Part: ALBR20SD4. HCPCS: E2615
Deep Contour, 16" Wide x 21" or 22" Overall Height..... \$735 Part: ALBR16TD4. HCPCS: E2615	Deep Contour, 20" Wide x 21" or 22" Overall Height..... \$735 Part: ALBR20TD4. HCPCS: E2615
Deep Contour, 16" Wide x 23" or 24" Overall Height..... \$735 Part: ALBR16XD4. HCPCS: E2615	Deep Contour, 20" Wide x 23" or 24" Overall Height \$735 Part: ALBR20XD4. HCPCS: E2615

1e. ADI Back, AL Series, 4-Point, Extra Deep Contour

Available with Back Cane Option only. Select a back type and the specific overall back height in subsequent section. Previously selected back cane height must be equal to the greater of overall height range.

Extra Deep Contour, 16" Wide x 21" or 22" Overall Height \$822 Part: ALBR16T64. HCPCS: E2620	Extra Deep Contour, 18" Wide x 23" or 24" Overall Height..... \$822 Part: ALBR18X64. HCPCS: E2620
Extra Deep Contour, 16" Wide x 23" or 24" Overall Height..... \$822 Part: ALBR16X64. HCPCS: E2620	Extra Deep Contour, 20" Wide x 19" or 20" Overall Height..... \$822 Part: ALBR20S64. HCPCS: E2620
Extra Deep Contour, 18" Wide x 19" or 20" Overall Height..... \$822 Part: ALBR18S64. HCPCS: E2620	Extra Deep Contour, 20" Wide x 21" or 22" Overall Height..... \$822 Part: ALBR20T64. HCPCS: E2620
Extra Deep Contour, 18" Wide x 21" or 22" Overall Height..... \$822 Part: ALBR18T64. HCPCS: E2620	Extra Deep Contour, 20" Wide x 23" or 24" Overall Height..... \$822 Part: ALBR20X64. HCPCS: E2620

1f. Q-Back 2.0, AL Series, Active Contour

Not available with Back Cane Option. Select a back type and the specific overall back height in subsequent section.

Active Contour, 14"Wx19-20"H Tilt/22"-23" H Recline..... \$568 Part: QB2AL14AS. HCPCS: E2613	Active Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$568 Part: QB2AL18AT. HCPCS: E2613
Active Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$568 Part: QB2AL14AT. HCPCS: E2613	Active Contour, 18"Wx23-24"H Tilt/26"-27" H Recline..... \$568 Part: QB2AL18AX. HCPCS: E2613
Active Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$568 Part: QB2AL16AS. HCPCS: E2613	Active Contour, 18"Wx25-26"H Tilt/28"-29" H Recline..... \$703 Part: QB2AL18AXX. HCPCS: E2613
Active Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$568 Part: QB2AL16AT. HCPCS: E2613	Active Contour, 20"Wx19-20"H Tilt/22"-23" H Recline..... \$568 Part: QB2AL20AS. HCPCS: E2613
Active Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$568 Part: QB2AL16AX. HCPCS: E2613	Active Contour, 20"Wx21-22"H Tilt/24"-25" H Recline..... \$568 Part: QB2AL20AT. HCPCS: E2613
Active Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL16AXX. HCPCS: E2613	Active Contour, 20"Wx23-24"H Tilt/26"-27" H Recline..... \$568 Part: QB2AL20AX. HCPCS: E2613
Active Contour, 18"Wx19-20"H Tilt/22"-23" H Recline..... \$568 Part: QB2AL18AS. HCPCS: E2613	Active Contour, 20"Wx25-26"H Tilt/28"-29" H Recline..... \$703 Part: QB2AL20AXX. HCPCS: E2613

1g. Q-Back 2.0, AL Series, Deep Contour

Not available with Back Cane Option. Select a back type and the specific overall back height in subsequent section.

Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$643 Part: QB2AL14DS. HCPCS: E2620	Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$643 Part: QB2AL18DT. HCPCS: E2620
Deep Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$643 Part: QB2AL14DT. HCPCS: E2620	Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$643 Part: QB2AL18DX. HCPCS: E2620
Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline \$643 Part: QB2AL16DS. HCPCS: E2620	Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL18DXX. HCPCS: E2620
Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline \$643 Part: QB2AL16DT. HCPCS: E2620	Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$643 Part: QB2AL20DS. HCPCS: E2620
Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$643 Part: QB2AL16DX. HCPCS: E2620	Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$643 Part: QB2AL20DT. HCPCS: E2620
Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL16DXX. HCPCS: E2620	Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$643 Part: QB2AL20DX. HCPCS: E2620
Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$643 Part: QB2AL18DS. HCPCS: E2620	Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$703 Part: QB2AL20DXX. HCPCS: E2620

TRU Balance 3 seating section continued on following pages.

10. TRU BALANCE 3 SEATING CONTINUED

1h. Q-Back 2.0, AL Series, Extra Deep Contour

Not available with Back Cane Option. Select a back type and the specific overall back height in subsequent section.

- | | |
|---|---|
| Extra Deep Contour, 14"Wx19-20"H Tilt/22"-23" H Recline \$703
Part: QB2AL14ES. HCPCS: E2620 | Extra Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline \$703
Part: QB2AL18ET. HCPCS: E2620 |
| Extra Deep Contour, 14"Wx21-22"H Tilt/24"-25" H Recline \$703
Part: QB2AL14ET. HCPCS: E2620 | Extra Deep Contour, 18"Wx23-24"H Tilt/26"-27" H Recline \$703
Part: QB2AL18EX. HCPCS: E2620 |
| Extra Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline..... \$703
Part: QB2AL16ES. HCPCS: E2620 | Extra Deep Contour, 18"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL18EXX. HCPCS: E2620 |
| Extra Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline..... \$703
Part: QB2AL16ET. HCPCS: E2620 | Extra Deep Contour, 20"Wx19-20"H Tilt/22"-23" H Recline \$703
Part: QB2AL20ES. HCPCS: E2620 |
| Extra Deep Contour, 16"Wx23-24"H Tilt/26"-27" H Recline \$703
Part: QB2AL16EX. HCPCS: E2620 | Extra Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$703
Part: QB2AL20ET. HCPCS: E2620 |
| Extra Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL16EXX. HCPCS: E2620 | Extra Deep Contour, 20"Wx23-24"H Tilt/26"-27" H Recline \$703
Part: QB2AL20EX. HCPCS: E2620 |
| Extra Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$703
Part: QB2AL18ES. HCPCS: E2620 | Extra Deep Contour, 20"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AL20EXX. HCPCS: E2620 |

1i. Q-Back 2.0, AFT Series, Active Contour

Not available with Back Cane Option. Select a back type and the specific overall back height in subsequent section.

- | | |
|--|--|
| AFT, Active, 16"Wx19-20"H Tilt/22"-23" H Recline \$568
Part: QB2AFT16AS. HCPCS: E2613 | AFT, Active, 18"Wx19-20"H Tilt/22"-23" H Recline \$568
Part: QB2AFT18AS. HCPCS: E2613 |
|--|--|

1j. Q-Back 2.0, AFT Series, Deep Contour

Not available with Back Cane Option. Select a back type and the specific overall back height in subsequent section.

- | | |
|--|--|
| AFT, Deep Contour, 16"Wx19-20"H Tilt/22"-23" H Recline..... \$643
Part: QB2AFT16DS. HCPCS: E2620 | AFT, Deep Contour, 18"Wx21-22"H Tilt/24"-25" H Recline..... \$643
Part: QB2AFT18DT. HCPCS: E2620 |
| AFT, Deep Contour, 16"Wx21-22"H Tilt/24"-25" H Recline..... \$643
Part: QB2AFT16DT. HCPCS: E2620 | AFT, Deep Contour, 20"Wx21-22"H Tilt/24"-25" H Recline \$643
Part: QB2AFT20DT. HCPCS: E2620 |
| AFT, Deep Contour, 18"Wx19-20"H Tilt/22"-23" H Recline \$643
Part: QB2AFT18DS. HCPCS: E2620 | |

1k. Q-Back 2.0, AFT Series, Extra Deep Contour

Not available with Back Cane Option. Select a back type and the specific overall back height in subsequent section.

- | |
|---|
| AFT, Extra Deep Contour, 16"Wx25-26"H Tilt/28"-29" H Recline \$703
Part: QB2AFT16EXX. HCPCS: E2620 |
|---|

10.12. Lateral Positioning Wedge

Optional.

- Lateral Positioning Wedge, pair \$16 ea.
Part: SET123465 x2
Available with Sport and TRU-Comfort backs. Not available on power tilt backs 16"-17"H with limited adjustability on 18" H. Not available on power recline back 19"-20"H with limited adjustability on 21"H.

10.13. Back Width Package

REQUIRED except when Omit Back Pan/Cushion is selected. Available with +/-4" of seat width.

- | | | |
|--|--------------------------------------|--------------------------------------|
| Match Seat Width Selections
Part: PTOINDV3468 | 14"W-16"W Package
Part: XRF000037 | 18"W-22"W Package
Part: XRF000039 |
| 10"W-14"W Package
Part: XRF000036 | 16"W-20"W Package
Part: XRF000038 | 20"W-24"W Package
Part: XRF000040 |

10.14. Back Width

REQUIRED except when Omit Back Pan/Cushion is selected. Part number is PTOINDV plus number following # below. Available with +/-4" of seat width.

- | | | | | | | | | | |
|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|
| 12" W
(#3562) | 14" W
(#3563) | 15" W
(#3856) | 16" W
(#3564) | 17" W
(#3565) | 18" W
(#3566) | 19" W
(#3567) | 20" W
(#3568) | 21" W
(#3569) | 22" W
(#3570) |
|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|

10.15. Overall Back Height

Static Seat and Power Tilt available height range: 14"-25"H. **Odd heights include 4" gap and even heights include 3" gap** between seat pan and bottom of back cushion.
Power Recline available height range: 19"-28"H. **Odd heights include 6" gap and even heights include 7" gap.** See the Sport and TRU-Comfort back availability chart for stock order back size ranges.

- | | | | | |
|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| 14" Height
Part: SETHGHT1027 | 17" Height
Part: SETHGHT1030 | 20" Height
Part: SETHGHT1033 | 23" Height
Part: SETHGHT1036 | 26" Height
Part: SETHGHT1039 |
| 15" Height
Part: SETHGHT1028 | 18" Height
Part: SETHGHT1031 | 21" Height
Part: SETHGHT1034 | 24" Height
Part: SETHGHT1037 | 27" Height
Part: SETHGHT1040 |
| 16" Height
Part: SETHGHT1029 | 19" Height
Part: SETHGHT1032 | 22" Height
Part: SETHGHT1035 | 25" Height
Part: SETHGHT1038 | 28" Height
Part: SETHGHT1041 |

TRU Balance 3 seating section continued on following pages.

10. TRU BALANCE 3 SEATING CONTINUED

Overall Height	Tilt or Static Width											Recline Width										
	12"	14"	15"	16"	17"	18"	19"	20"	21"	22"	12"	14"	15"	16"	17"	18"	19"	20"	21"	22"		
	14"	T	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
15"	T	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-		
16"	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T		
17"	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T		
18"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	
19"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	
20"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	
21"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	
22"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	
23"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	
24"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	
25"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	
26"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	
27"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	
28"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	
29"	-	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	S	T	
30"	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
31"	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
32"	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	

Sport Back Q-Back TRU-Comfort Back

Sport, Q-Back, and TRU-Comfort Back Availability Chart

10.16. Power Tilt Privacy Flap

Optional. A privacy flap is standard with power recline systems. This selection is to add a flap for power tilt only systems.

Privacy Flap \$48
 Part: XRF000351

10.17. Back Height Extension

Optional. Not available with back cane option, 12" back width, or 24" back width. Adds 2" of seat back pan height. 25"H selection includes a back height extension.

Back Height Extension \$48
 Part: XRF000041

11. TRU BALANCE 3 ARMRESTS

REQUIRED. Select the desired type of armrest below. Please read notes to ensure seating compatibility. Be sure to also select an armrest size and armpad type.

11.1a. Flip-Back, 8"-16" H Armrests

Choose a left and right armrest and armrest height in the following section. Not available with waterfall armpads. Not available with Back Cane Option.

Flip-Back, 8"-16"H Armrests Standard Flip-Back, 8"-16"H Armrests for Static Seat \$289
 Left part: ARM1806662. Right part: ARM1806663 Left part: XRF000320. Right part: XRF000321
 Not available with TRU Balance 3 Static Seat.

Armrest Support

Optional. Available up to 2" difference between seat and back width via endomorph or mesomorph configuration. Not available less than 16" width seat or back.

Flip-Back Armrest Downpost \$92 ea.
 Left part: XRF000265. Right part: XRF000266

11.1b. Drop-In, Height Adjustable, Removable, T-Arm Armrests

Only available with TRU Balance 3 Static Seat. Includes flat armpads.

Full Length, Drop-In Armrests No Charge Desk Length, Drop-In Armrests No Charge
 Left part: ARM1803253. Right part: ARM156742 Left part: ARM1803254. Right part: ARM156743

11.1c. Two-Post, Height Adjustable, Removable Armrests

Select the armrest and the left and right armrest heights in the following section. These armrests may interfere with offset backs/ back canes.

1a. Two-Post Armrest for Power Seating

Full Length, Two-Post Armrest for Power Seating No Charge Desk Length, Two-Post Armrest for Power Seating.... No Charge
 Left part: XRF000047. Right part: XRF000048 Left part: XRF000049. Right part: XRF000050

1b. Two-Post Armrest for Static Seating

Full Length, Two-Post Armrest for Static Seating.....\$158 ea. Desk Length, Two-Post Armrest for Static Seating\$158 ea.
 Left part: XRF000051. Right part: XRF000052. HCPCS: E0973 Left part: XRF000053. Right part: XRF000054. HCPCS: E0973

11.1d. Heavy Duty, Removable Armrests

HD removable armrests require armrest height and armrest size selections. Only one height package is selectable as a normal configuration. Payer coding, coverage, and bundling guidelines may apply.

8"-12"H, Single Post, Height Adjustable, Removable, \$158 Heavy Duty (Sold as pair)
 Part: XRF000055. HCPCS w/ Static Seat: E0973

10"-14"H, Single Post, Height Adjustable, Removable,..... \$158 Heavy Duty (Sold as pair)
 Part: XRF000056. HCPCS w/ Static Seat: E0973

TRU Balance 3 armrests section continued on following page.

11. TRU BALANCE 3 ARMRESTS CONTINUED

11.1e. Pediatric, Removable, Height Adjustable Armrests

Height adjustment from 7.75" to 9.75" in .5" increments. **Not available with full length armpads.**

1. Pediatric Armrest Selection

L R Pediatric, Removable Armrests \$79 ea.
 Left Part: ARM142876. Right part: ARM142877. HCPCS: E0973

2. Pediatric, Left Armrest Height

7.75" H 8.25" H
 PTOINDV3722 PTOINDV3724

3. Pediatric, Right Armrest Height

7.75" H 8.25" H
 PTOINDV3721 PTOINDV3723

11.2. Armrest Height

Do not select with Pediatric Armrests. Set left and right armrest height. Note the available height range in the option description.

1. Left Armrest Height

8" H PTOINDV3535	10" H PTOINDV3537	12" H PTOINDV3539	14" H PTOINDV3541	16" H PTOINDV3543
9" H PTOINDV3536	11" H PTOINDV3538	13" H PTOINDV3540	15" H PTOINDV3542	

2. Right Armrest Height

8" H PTOINDV3450	10" H PTOINDV3452	12" H PTOINDV3454	14" H PTOINDV3456	16" H PTOINDV3458
9" H PTOINDV3451	11" H PTOINDV3453	13" H PTOINDV3455	15" H PTOINDV3457	

11.3. Armrest Size

Select a left (L) and right (R) armrest size. This also affects armpad selection.

<input type="checkbox"/> L <input type="checkbox"/> R 10" Length Left part: XRF000057. Right part: XRF000123	<input type="checkbox"/> L <input type="checkbox"/> R 14" Length Left part: XRF000058. Right part: XRF000124	<input type="checkbox"/> L <input type="checkbox"/> R 17" Length Left part: XRF000059. Right part: XRF000125 Only available with Flip-Back, 8"-16" H Armrests and straight arm pads.
--	--	--

12. TRU BALANCE 3 ARMPADS

REQUIRED. Select to omit an armpad selection or select a left and right armpad type. Please note that waterfall and gel armpads may cause interference with the seat back.

12.1a. Omit Armpad

Select one below if aftermarket armpads or handpads are required. See [Positioning Components Order Form](#).

No Armpad Selection..... No Charge Part: XRF000060 Not available with Flip-back, 8"-16" H armrests.	Interface for Aftermarket Pad No Charge Part: XRF000308 Only available with Flip-back, 8"-16" armrests.
---	---

12.1b. Foam Armpads

<input type="checkbox"/> L <input type="checkbox"/> R Straight, Armpad Standard Left part: XRF000061. Right part: XRF000126	<input type="checkbox"/> L <input type="checkbox"/> R Waterfall, Armpad No Charge Left part: XRF000062. Right part: XRF000127 Not available with Flip-back, 8"-16" armrests.
--	--

12.1c. Gel Armpads

Sold only in pairs. Armpads are 3" wide. Manufacturer part number for full length is 314G and desk length is 310G.

a. Falcon Gel Armpads for Flip-Back Armrests

Falcon Gel, Full Length Left and Right\$195 Part: ARM1806807	Falcon Gel, Desk Length Left, Full Length Right.....\$171 Part: ARM1806810
Falcon Gel, Full Length Left, Desk Length Right.....\$171 Part: ARM1806813	Falcon Gel, Desk Length Left and Right\$147 Part: ARM1806806

b. Falcon Gel Armpads for Tube Style Armrests

Falcon Gel, Full Length Left and Right\$311.22 Part: ARM136438	Falcon Gel, Desk Length Left, Full Length Right.....\$287.93 Part: ARM136439
Falcon Gel, Full Length Left, Desk Length Right.....\$287.93 Part: ARM136442	Falcon Gel, Desk Length Left and Right \$265.02 Part: ARM136437

13. TRU BALANCE 3 JOYSTICK BRACKETS

REQUIRED.

Standard Fixed Receiver Standard Part: XRF000063	Stealth Height Adjustable Receiver..... \$236 Part: XRF000064. HCPCS: K0108
---	--

14. TRU BALANCE 3 HEADRESTS

Optional. Select a headrest pad and bracket.

14.1. Headrest, Comfort Plus with Cool Core

1. Headrest Pads

Stealth 8" Comfort Plus, Cool Core (CP08CCFPM)\$179
 Part: ST-CP08CCFPM-Q. HCPCS: E0955

Stealth 14" Comfort Plus, Cool Core (CP14CCFPM)\$179
 Part: ST-CP14CCFPM-Q. HCPCS: E0955

Stealth 10" Comfort Plus, Cool Core (CP10CCFPM).....\$179
 Part: ST-CP10CCFPM-Q. HCPCS: E0955

Stealth 6", Comfort Plus, Cool Core (CP06CCFPM)\$179
 Part: ST-CP06CCFPM-Q. HCPCS: E0955

2. Headrest Brackets

Comfort Plus Headrest Mounting Bracket, UniLink..... \$311
 Part: ST-UHB2D33200-Q. HCPCS: E1028

15. LAP BELTS

REQUIRED. Select a lap belt below.

15.1a. Standard Lap Belts

38" Lap Belt Part: ACCBELT1006	50" Lap Belt Part: ACCBELT1000	60" Lap Belt Part: ACCBELT1002	70" Lap Belt Part: ACCBELT1003	80" Lap Belt Part: ACCBELT1004	90" Lap Belt Part: ACCBELT1001
--	--	--	--	--	--

15.1b. Optional Lap Belts

Stealth Retractable Lap Belt (IFRB5350).....\$181
 Part: ST-IFRB5350

Stealth Positioning 60" Single Pull, Push Button, Extra Large.....\$50
 Padded (11.5"), 2" Belt (IFPB5315)
 Part: ST-IFPB5315-Q

Stealth Positioning 60" Single Pull, Push Button, Medium..... \$50
 (7.5") Padded, 1.5" Belt (IFPB5213)
 Part: ST-IFPB5213-Q

60" Sure-Lok Lap Belt (Non-retractable)..... \$120
 Required with ISO 7176-19 Compliant Transit Option.
Coming soon. Part: ACCBELT1016

Stealth Positioning 60" Single Pull, Push Button, Large \$50
 Padded (9.5"), 2" Belt (IFPB5314)
 Part: ST-IFPB5314-Q

16. POWER LEGREST CONTROLS

Optional. Specific power legrest options are selected later in this form. Select the control type here, if power legrests are desired, to assist the selection of the power positioning electronics.

Independent Legrests (Requires 2 actuator controls)
 Part: XRF000108

Articulating Foot Platform (Requires 1 Actuator Control)
 Part: XRF000110

Combined Legrests (Requires 2 actuator control)
 Part: XRF000109

17. LEGRESTS

Optional. All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit.

17.1a. Basic Center Mount Foot Platform

Select legrest length in following section.

Short	Medium	Long
8"-10"	10"-14"	14"-19"
Basic Center Mount Extension Range		

1. Center Mount Foot Platform Type

Basic Center Mount..... No Charge
 Part: XRF000323

2. Basic Center Mount Foot Platform Extension and Angle Adjustable Footplates

See legrest length matrix for compatibility.

Short Extension, 5"x9" FootplatesNo Charge
 Part: MEC1906675

Medium Lower Extension, 6.25"W Tapered 5"W x 12"D... No Charge
 Part: MEC2104853

Medium Extension, 5"x9" FootplatesNo Charge
 Part: MEC1906680

Long Lower Extension, 6.25"W Tapered 5"W x 12"D..... No Charge
 Part: MEC2104854

Long Extension, 5"x9" FootplatesNo Charge
 Part: MEC1906677

3. Basic Center Mount Calf Pads

Small Calf Pads, Basic Center Mount..... \$158
 Part: RIG2104369

Large Calf Pads, Basic Center Mount \$158
 Part: RIG2202997
 Only available with Long Length Extension.

4. Basic Center Mount Additional Accessories

Basic Center Mount Footplate Connector \$16
 Part: RIG2104240

17.1b. Center Mount Foot Platform

If power articulating foot platform is selected, be sure power legrest control has been selected in previous section. Select legrest length in following section.

1. Center Mount Foot Platform Type

Power Articulating Foot Platform \$3,145
 Part: PTOINDV3754. HCPCS: E1012+K0040

Legrests section continued on following pages.

17. LEGRESTS CONTINUED

2. Extension and Angle Adjustable Footplates

See legrest length matrix for compatibility.

Medium Extension, 5"x6" Footplates
Part: RIG150206

Medium Extension, 5"x8" Footplates
Part: RIG150203

Medium Extension, 5"x11" Footplates
Part: RIG150212

Long Extension, 5"x6" Footplates
Part: RIG151505

Long Extension, 5"x8" Footplates
Part: RIG151506

Long Extension, 5"x11" Footplates
Part: RIG151508

Medium	Long
11"-17"	14"-20"
Power AFP Extension Range	

3a. Center Mount Calf Support

*Calf pads are standard with a power articulating foot platform. Select calf pads covers for additional charge. Calf pad and panel covers include 1/2" of soft visco foam.

Small Calf Pads, AFP Only, 6"Wx6"H (pair) Standard
Part: RIGASMB7120018

Large Calf Pads, AFP Only, 8" taper to 6"Wx8"H (pair) Standard
Part: RIGASMB7120019

3b. Center Mount Calf Panel

Calf pads are standard with a power articulating foot platform.

Small Calf Panel for AFP, 18"Wx8"H..... \$184
Part: FRMASMB7110007

Large Calf Panel for AFP, 24"Wx8"H \$237
Part: RIG123824

Only available on widths of 20" or greater.

4. Modular Footplate Extension Components

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

Width Extension for 8"D Footplates \$69
Part: RIG152517

Width Extension for 11"D Footplates \$69
Part: RIG152518

2" Bell-Shaped Length Extension \$69
Part: RIG152519
Not compatible with width extensions. Can be used with 6" outside edge with 11" footplates.

4" Outside Edge (pair) \$69
Part: RIG152510

Only available with 6"D footplates and 8"D footplate extension.

6" Outside Edge (pair) \$69
Part: RIG152516
Available with 8"D and 11"D footplates and 11"D footplate extension.

5. Additional Accessories

Small Reversed Dartex Calf Pad Cover (pair) \$158
Part: RIG133717

Large Reversed Dartex Calf Pad Cover (pair) \$158
Part: RIG133719

Small Reversed Dartex Panel Cover (pair)..... \$210
Part: RIG134673

Large Reversed Dartex Panel Cover (pair) \$210
Part: RIG134744

Small Foam Calf Pad Cover (pair)..... \$67
Part: ST-CSCS0775-Q

Large Foam Calf Pad Cover (pair) \$83
Part: ST-CSCS0985-Q

Small Foam Panel Cover (pair) \$176
Part: ST-CSCP0920-Q

Large Foam Panel Cover (pair)..... \$282
Part: ST-CSCP0927-Q

Footplate Connector \$16
Part: RIGASMB7120023

Footplate Riser \$75
Part: XRF000348

Not available with Long Extension. Select legrest length between 2" to 9" at 1" increments.

Footplate Riser Calf Panel \$50
Part: FRMASMB7110011
Only available and required with Footplate Riser selection and extension lengths of 6" to 9".

6. Stealth Gel Footplate Covers

5"x8" Tapered Gel Footplate Covers \$222
Part: ST-FPG5008T-Q. HCPCS: K0108

5"x11" Tapered Gel Footplate Covers..... \$243
Part: ST-FPG5011T-Q. HCPCS: K0108

5"x8" Tapered Gel Footplate Covers for Width Extension \$243
Part: ST-FPG6508T-Q. HCPCS: K0108

5"x11" Tapered Gel Footplate Cover for Width Extension \$270
Part: ST-FPG7011T-Q. HCPCS: K0108

17.1c. Power Elevating Legrests

Be sure power legrest control has been selected in previous section. Select an upper extension type and footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

1. Power Elevating Legrest Upper Assembly

Includes calf pads.

Power Elevating Legrests, Pair \$2,516
Part: XRF000341. HCPCS: E1010+K0040+E0951

Legrest Length Matrix	Seat to Floor Option		
	Low	Middle	High
Power ELR: Static Seat, Recline Only			
8.5"-14.5" Lower Length	8.5"-12.5"	8.5"-13.5"	8.5"-14.5"
12.5"-21" Lower Length	N/A	12.5"-13.5"	12.5"-14.5"

Legrest Length Matrix	Seat to Floor Option		
	Low	Middle	High
Power ELR: Tilt Only, Tilt and Recline, Any Power Lift			
8.5"-14.5" Lower Length	8.5"-14.5"	8.5"-14.5"	8.5"-14.5"
12.5"-21" Lower Length	12.5"-15.5"	12.5"-15.5"	12.5"-17"

2. Lower Extension

8.5" to 14.5" Lower Length
Part: XRF000344

12.5" to 21" Lower Length
Part: XRF000347

Legrests section continued on following pages.

17. LEGRESTS CONTINUED

3. Angle Adjustable Footplate with Heel Loop and Wheel Bumper

*The 6"W footplates may cause interference with seat widths less than 16". *The 8"W footplates may cause interference with seat widths less than 19".

4"Wx6"D, Pair Part: FRMASMB7448	5"Wx8"D, Pair Part: FRMASMB7450	6"Wx11"D ¹ , Pair Part: FRMASMB7452
5"Wx6"D, Pair Part: FRMASMB7449	6"Wx8"D ¹ , Pair Part: FRMASMB7451	8"Wx11"D ² , Pair Part: FRMASMB7453

17.1d. Heavy Duty, Drop-In Legrests

Select a configuration, a legrest, extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

1. Drop-In Legrest Selection

<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R	70° HD Drop-in for 8.5"-12" Lower Extension No Charge Pair part: FRMASMB4762. Left: RIG137706. Right: RIG137707	<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R	70° HD Drop-in 13"-18" or 16"-21" Lower No Charge Extension Pair part: FRMASMB4853. Left: RIG137708. Right: RIG137709
--	--	--	---

2. Lower Extension

<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R	8.5"-12" HD Drop-Ins/12.5"-16.5" ELR, Ext Length Pair part: XRF000101. Left: XRF000102. Right: XRF000103.	<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R	16"-21" HD Drop-Ins/18.5"-22.5" ELR, Ext Length Pair part: FRMASMB4754. L: FRMASMB4725. R: FRMASMB5393
<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R	13"-18" HD Drop-Ins/15.5"-19.5" ELR, Ext Length Pair part: XRF000104. Left: XRF000105. Right: XRF000106		

3. Angle Adjustable Footplate(s) with Heel Loops

<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R	6"Wx8"D ¹ \$90 ea. Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795. HCPCS: K0040+E0951	<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R	8"Wx11"D ² \$90 ea. Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799. HCPCS: K0040+E0951 May cause interference on seat widths less than 20".
<input type="checkbox"/> P <input type="checkbox"/> L <input type="checkbox"/> R	6"Wx11"D ¹ \$90 ea. Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797. HCPCS: K0040+E0951 May cause interference on seat widths less than 18".		

17.2. Set Legrest Extension Length

REQUIRED with matching left and right legrest configuration. Part number is PTOINDV plus number following # below. Non-matching or individual legrest configurations lengths are not set by the manufacturer.

1a. Centermount Extension Length

2" L (#3989)	5" L (#3992)	8" L (#3680)	9.5" L (#3683)	11" L (#3686)	12.5" L (#3689)	14" L (#3692)	15.5" L (#3695)	17" L (#3698)	18.5" L (#3701)
3" L (#3990)	6" L (#3993)	8.5" L (#3681)	10" L (#3684)	11.5" L (#3687)	13" L (#3690)	14.5" L (#3693)	16" L (#3696)	17.5" L (#3699)	19" L (#3702)
4" L (#3991)	7" L (#3994)	9" L (#3682)	10.5" L (#3685)	12" L (#3688)	13.5" L (#3691)	15" L (#3694)	16.5" L (#3697)	18" L (#3700)	

1b. Front Rigging Extension Length

Non-matching or individual legrest configurations lengths are not set by the manufacturer. See previous section for center mount extension lengths.

7" L (#3470)	8.5" L (#3473)	10" L (#3476)	11.5" L (#3479)	13" L (#3482)	14.5" L (#3485)	16" L (#3488)	17.5" L (#3491)	19" L (#3494)	20.5" L (#3497)	22" L (#3669)
7.5" L (#3471)	9" L (#3474)	10.5" L (#3477)	12" L (#3480)	13.5" L (#3483)	15" L (#3486)	16.5" L (#3489)	18" L (#3492)	19.5" L (#3495)	21" L (#3498)	22.5" L (#3670)
8" L (#3472)	9.5" L (#3475)	11" L (#3478)	12.5" L (#3481)	14" L (#3484)	15.5" L (#3487)	17" L (#3490)	18.5" L (#3493)	20" L (#3496)	21.5" L (#3600)	

18. TRU BALANCE 3 POWER POSITIONING ELECTRONICS

REQUIRED with power seating or legrest options.

18.1a. Actuator Control through Joystick

1. Actuator Control Selection

Single Actuator Through Controller.....\$2,163 Interface Part: XRF000335. HCPCS: E2310	Multiple Actuator Through Controller and \$3,040 Interface Part: XRF000336. HCPCS: E2311
--	--

2. Q-Logic 3 Optional Switch Control

Only available with Q-Logic 3, multiple actuator configurations.

4-Function Push Button Part: CTL123218 See section 21.2 for switch mounting.	4-Way Toggle Part: ELEASMB5925	4-Way Button Toggle Part: ELEASMB5926
--	-----------------------------------	--

3. User Controls Side Selection

Required with 4-Way Toggle and 4-Way Button Toggle.

Mount User Controls on Left Part: PTOINDV3996	Mount User Controls on Right Part: PTOINDV3997
--	---



TRU Balance 3 power positioning electronics section continued on following page.

18. TRU BALANCE 3 POWER POSITIONING ELECTRONICS CONTINUED

18.1b. Actuator Control with iAccess Switch

Single Actuator Control Through Flat Button.....\$2,547 iAccess and Controller and Interface Part: XRF000337. HCPCS: E2310	Multiple Actuator Control Through Flat Button\$3,460 iAccess and Controller and Interface Part: XRF000339. HCPCS: E2311
Single Actuator Control Through Toggle Button\$2,547 iAccess and Controller and Interface Part: XRF000338. HCPCS: E2310	Multiple Actuator Control Through Toggle Button.....\$3,460 iAccess and Controller and Interface Part: XRF000340. HCPCS: E2311

18.2. 4-Function Push Button/iAccess Mount

REQUIRED with 4-Function Push Button or iAccess options. Under arm mount is not available with TRU Balance 3 HD armrests.

Joystick Mount No Charge Part: MEC1907500	Attendant Mount \$26 Part: MEC1906564
Under Arm Mount, Flip-Back Armrest..... No Charge Part: MEC1906531	Under Arm Mount, Flip-Back Armrest, plus iAccess Guard \$26 Part: MEC2001400
Under Arm Mount, Tube Style Armrest No Charge Part: MEC1906530	Under Arm Mount, Tube Style Armrest, plus iAccess Guard \$26 Part: MEC2001401

1. Under Arm Mount Side Selection

Left Side Under Arm Mount Part: PTOINDV3678	Right Side Under Arm Mount Part: PTOINDV3679
--	---

19. TRU BALANCE 3 ATTENDANT CONTROL


Optional. *Attendant control will be at no charge when ordered with Tilt and Recline, alternative drive controls, and a Q-Logic Display or SCIM module. Mounting bracket is included.

Attendant Control \$525 Part: XRF000078. HCPCS: E2331 Select this option if no charge option criteria is not met.	Attendant Control with Tilt, Recline, Alternative No Charge Controls and Q-Logic Display or SCIM Module Part: XRF000079
---	---

20. PROVIDER INSTALLED ACCESSORIES

Optional. Multiple accessories can be mounted along the various tracks but some interference may occur.

TRU Balance 3 Cane/Crutch Holder \$105 Part: ACC152515. HCPCS: E2207	TRU Balance 3 Oxygen Holder \$210 Part: ACC125006. HCPCS: E2208 Includes accessory bar, a credit will be applied PTOIN-DV3769.
Quantum Backpack \$168 Part: ACC1905853	TRU Balance 3 Clothing Guard (Pair) \$237 Part: ACC125007 Not available with Heavy Duty armrests.
Grip Lap Board (20005)..... \$38.95 Part: ACC152992	TRU Balance 3 Glove Box (C) \$158 Part: ACC125011
Grip Drink Holder (20085) \$15.50 Part: ACC152993	Stealth Large Tablet Holder \$420 Part: ACC141818 Not available with Heavy Duty armrests.
TRU Balance 3 Hydration System (B) \$237 Part: ACC125010	Compact Vent Tray \$1,465 Part: ACC133889. HCPCS: E1029
TRU Balance 3 Transfer Handle (D) \$105 ea. Part: MEC1809930. HCPCS: E0705 Not available with Heavy Duty armrests.	Push Handles (pair)..... \$195 Part: ACC2005027 Not available with oxygen holder, back cane option, or vent tray.
TRU Balance 3 Cup Holder No Charge Part: ACC125003 For use with TRU Balance 3 (track style) armrests.	Accessory Bar \$132 Part: XRF000128 Horizontal rail attached to rear of seat. Standard with the Trilogy vent tray and oxygen holder. Not available with back cane option.
Cup Holder (Tube-style Armrest)..... No Charge Part: ACCASMB2634 For use with tube-style armrests.	Large, Black, Weather Cover \$202 Part: ACCCOVR1021
TRU Balance 3 Personal Item Hook \$69 Part: ACC125008 Not available with back cane option.	Q-Logic EX/NE Programmer \$573 Part: ELEASMB6976
TRU Balance 3 Reflector Kit \$53 Part: ACC124670	Flip-Back Armrest, Side Track Accessory Bracket ... \$11 Part: MEC1809178 Only available with flip-back, 8"-16" H armrests. This adds a 47mm channel to the side of the armrest to allow mounting of accessories with t-nuts.
TRU Balance 3 Backpack Holder (Pair) (A) \$69 Part: ACC147774 Not available with back cane option.	Privacy Bag \$42 Part: ACC1907570
TRU Balance 3 Front Mounted \$69 Cell Phone/ Media Holder, X-Grip (E) Part: ACC155835 Not available with Heavy Duty armrests.	Personal Item Pouch \$84 Part: ACC2101017
TRU Balance 3 Rear View Mirror..... \$69 Part: ACC155833 Not available with Heavy Duty armrests.	



21. BRANCH PATCHES

Optional.

Air Force Patch\$21 Part: ACC140240	Coast Guard Patch\$21 Part: ACC140239	Navy Patch\$21 Part: ACC140244
Army Patch\$21 Part: ACC140241	Marines Patch\$21 Part: ACC140243	Purple Heart Patch\$21 Part: ACC140238



22. POSITIONING COMPONENTS

Optional. More options are available on the [Positioning Components Order Form](#).

22.1. Stealth Thoracic Laterals

Swing-away laterals may interfere with TRU Balance 3 armrests and 2-post, flip up armrests.

1. Bracket

Swing-away, UniLink, Redesign Back, pair (ULS4042000) \$664 Part: ST-ULS4042000-Q. HCPCS: E1028 Only available with Redesigned back assemblies.	Swing-away, UniLink, Q-Back 2, Left (ULPK404103L00) \$332 Part: ST-ULPK404103L00-Q. HCPCS: E1028 Only available with Active Contour Q-Backs.
Swing-away, UniLink, 125° w/Track, pair (ULPK404108B00) ... \$664 Part: ST-ULPK404108B00-Q. HCPCS: E1028 Only available with Sport or TRU Comfort backs.	Swing-away, UniLink, Q-Back 2, Right (ULPK404103R00) \$332 Part: ST-ULPK404103R00-Q. HCPCS: E1028 Only available with Active Contour Q-Backs.
Swing-away, UniLink, Back Cane, pair (ULPK404104B00) \$664 Part: ST-ULPK404104B00-Q. HCPCS: E1028 Only available with back cane and omit back configurations.	

2. Pads

Curved Pads, 3"x5", pair (TWBL-3X5)..... \$250 Part: INDPART2663. HCPCS: E0956	Curved Pads, 4"x6", pair (TWBL-4X6)..... \$250 Part: INDPART2664. HCPCS: E0956
---	---

3. Installation

Install Laterals Part: PTOINDV3499	Do Not Install Part: PTOINDV3500
---------------------------------------	-------------------------------------

22.2. Stealth Pelvic/Thigh Guides

Pelvic/thigh guides may interfere with the joystick receiver block and down posts. Please provide quantity next to check box.

1. Pelvic/Thigh Guide Bracket

Flip Down Bracket (TWBTLADD-HW) \$281 ea. Part: INDPART2804. HCPCS: E1028	Removable Bracket with 1" Offset (LP430 & LHW-121) \$382 ea. Part: INDPART2854. HCPCS: E1028
Removable Bracket (LP430)..... \$284 ea. Part: INDPART2805. HCPCS: E1028	

2. Pelvic/Thigh Guide Adapter

2" Height Adapter Plate (LHW-122)..... \$28 Part: INDPART2821. HCPCS: K0108 Height adapter is needed with TRU-Comfort seat cushions. Not available with the flip down bracket.	2 1/4" Adapter Plate (LHW-129) \$57 Part: POS140430. HCPCS: K0108 Only available with the flip down full surface bracket.
--	---

3. Pelvic/Thigh Guide Pad

^Please verify the appropriate billing code to be used for this item with all third-party payors.

Short Pad, 4"x4" (TWBADD-4X4)..... \$81 ea. Part: INDPART2801. HCPCS: E0953/E0956^	Long Pad, 4"x12" (TWBADD-4X12) \$81 ea. Part: INDPART2803. HCPCS: E0953/E0956^
Medium Pad, 4"x8" (TWBADD-4X8)..... \$81 ea. Part: INDPART2802. HCPCS: E0953/E0956^	

4. Pelvic/Thigh Guide Installation

Install on Left Part: PTOINDV3585	Install on Right Part: PTOINDV3586	Do Not Install Part: PTOINDV3719
--------------------------------------	---------------------------------------	-------------------------------------

22.3. Stealth UniLink Pelvic/Thigh Guide

1. UniLink Pelvic/Thigh Hardware

UniLink, 5" Height Adjustable & Removable Hardware \$265 ea. Part: ST-UTTP5-Q. HCPCS: E1028	UniLink, 7" Height Adjustable & Removable Hardware \$276 ea. Part: ST-UTTP7-Q. HCPCS: E1028
Flip Down with LEFT Full Surface Contact \$281 ea. Feature (TWBADD-L & LPHW407) Part: POS140428. HCPCS: E1028	Flip Down with RIGHT Full Surface Contact..... \$281 ea. Feature (TWBADD-R & LPHW407) Part: POS140429. HCPCS: E1028

2. UniLink Pelvic/Thigh Pads

^Please verify the appropriate billing code to be used for this item with all third-party payors.

UniLink, 4" x 4" Pad \$106 ea. Part: ST-UTP544056100-Q. HCPCS: E0953/E0956	UniLink, 4" x 12" Pad \$106 ea. Part: ST-UTP54D056100-Q. HCPCS: E0953/E0956
UniLink, 4" x 8" Pad..... \$106 ea. Part: ST-UTP548056100-Q. HCPCS: E0953/E0956	

Positioning components section continued on following page.

22. POSITIONING COMPONENTS CONTINUED

3. UniLink Pelvic/Thigh Additional Links

UniLink, 1.5" Link\$81 ea. Part: ST-UIL001150A000-Q QTY	UniLink, 2" Link.....\$81 ea. Part: ST-UIL001200A000-Q QTY
---	--



UniLink Pad & Hardware

4. UniLink Pelvic/Thigh Guide Installation

Install on Left Part: PTOINDV3585	Install on Right Part: PTOINDV3586	Do Not Install Part: PTOINDV3719
--------------------------------------	---------------------------------------	-------------------------------------

22.4. Stealth Positioning Chest Harness

Upper Classic, Medium (IFPS1103) \$154 ea. Part: ST-IFPS1103-Q. HCPCS: E0960	Upper Contour, Medium (IFPS2103)..... \$154 ea. Part: ST-IFPS2103-Q. HCPCS: E0960
Upper Classic, Large (IFPS1104) \$154 ea. Part: ST-IFPS1104-Q. HCPCS: E0960	Upper Contour, Large (IFPS2104) \$154 ea. Part: ST-IFPS2104-Q. HCPCS: E0960



Contoured and Classic Chest Harness

22.5. Stealth Positioning Chest Straps

Chest Strap, Medium w/ Belt Anchor\$97 ea. (IFCS4113) Part: ST-IFCS4113-Q. HCPCS: E0960	Chest Strap, Large w/ Belt Anchor\$97 ea. (IFCS4114) Part: ST-IFCS4114-Q. HCPCS: E0960
---	--



Chest Strap

22.6. Stealth Swing-Away Headrest Mounts

These swing-away headrest mounts are only available with the Stealth Comfort Series headrests.

1. Facial Lateral Hardware


Facial Lateral for Comfort Series..... \$318 ea. Right Side (CPS9-R-HW) Part: INDPART2647. HCPCS: E1028	UniLink Facial Lateral, Swing Away, Left \$318 Part: ST-SAR2A-L-Q. HCPCS: E1028
Facial Lateral for Comfort Series..... \$318 ea. Left Side (CPS9-L-HW) Part: INDPART2648. HCPCS: E1028	UniLink Facial Lateral, Dual Swing Away\$492 2 Upper & 1 Lower Right Part: ST-SAR4A-DBU-RL-Q. HCPCS: E1028
Facial Laterals for Comfort Series \$376 ea. Left and Right Side (CPS9-B-HW) Part: INDPART2649. HCPCS: E1028	UniLink Facial Lateral, Dual Swing Away\$492 2 Upper & 1 Lower Left Part: ST-SAR4A-DBU-LL-Q. HCPCS: E1028
UniLink Facial Lateral, Swing Away, Right..... \$318 Part: ST-SAR2A-R-Q. HCPCS: E1028	UniLink Facial Lateral, Dual Swing Away\$434 2 Upper Part: ST-SAR4A-DBU-Q. HCPCS: E1028



Left and Right Facial Lateral with Pads

2. Facial Lateral Pads and Accessories

Egg Switch Mount (ESM) (A)\$57 ea. Part: INDPART2815 QTY	Facial Pad, 6"x2", Cool Core (920M-Q).....\$61 ea. Part: ST-920M-Q. HCPCS: K0108 QTY
Figure 8 Bracket with 13" Rod..... \$289 ea. and Egg Switch Mount (B) Part: ST-F8675-13-Q QTY	Swing-Away Joystick Mount, 18" Rod (C) \$276 ea. Part: ST-SUS18-Q QTY Allows mounting of facial pads and attaches to the headrest
Facial Pad, 4".5x2", Cool Core (900M-Q)\$61 ea. Part: ST-900M-Q. HCPCS: K0108 QTY	




Headrest Accessories

22.7. Stealth Elbow Support

1. Elbow Support Bracket

Left, Swing-Away, Elbow Support \$199 ea. Bracket Tube Style Arms (ARMS-ES-L) Part: POS124064. HCPCS: E1028	Right, Swing-Away, Elbow Support..... \$199 ea. Bracket Tube Style Arms (ARMS-ES-R) Part: POS124066. HCPCS: E1028
Left, UniLink, Fixed, Elbow Support Bracket.....\$211 for Flip-Back, 8"-16"H Armrests (UETB4021FL-Q) Part: ST-UETB4021FL-Q	Right, UniLink, Fixed, Elbow Support Bracket.....\$211 for Flip-Back, 8"-16"H Armrests (UETB4021FR-Q) Part: ST-UETB4021FR-Q



Elbow Supports

2. Elbow Support Pad

3"Wx3"H Elbow Pad (ARMS-ES-3X3)\$57 ea. Part: POS124065. HCPCS: K0108 QTY	4"Wx4"H Elbow Pad (ARMS-ES-4X4).....\$57 ea. Part: POS124067. HCPCS: K0108 QTY
---	--

22.8. Stealth Foot Boxes

Stealth Small Foot Boxes, 5"Wx8"Dx7"H..... \$539 (SPLE107) (pair) Part: INDPART2786. HCPCS: E0954	Stealth Large Foot Boxes, 7"Wx12"Dx11"H..... \$539 (SPLE109) (pair) Part: INDPART2788. HCPCS: E0954
Stealth Medium Foot Boxes, 6"Wx10"Dx9"H \$539 (SPLE108) (pair) Part: INDPART2787. HCPCS: E0954	

22.9. Stealth Positioning Shoe Holder

Small Shoe Holder, 3/4" Webbing (IFSH4211).....\$87 ea. Part: ST-IFSH4211-Q. HCPCS: E0951+E0952 QTY	Large Shoe Holder, 3/4" Webbing (IFSH4213)\$87 ea. Part: ST-IFSH4213-Q. HCPCS: E0951+E0952 QTY
Medium Shoe Holder, 3/4" Webbing (IFSH4212)\$87 ea. Part: ST-IFSH4212-Q. HCPCS: E0951+E0952 QTY	

Positioning components section continued on following page.

23. LAP TRAY

23.1. Quantum Lap Tray

Mounting bracket is required with lap tray selection. These trays are compatible with Flip-Back armrests.

1a. Black Lap Tray

- Black Lap Tray, 11" W \$200
Part: MEC1905328-11B. HCPCS: E0950
For use with 12"-13" back width.
- Black Lap Tray, 13" W \$200
Part: MEC1905328-13B. HCPCS: E0950
For use with 14"-15" back width.
- Black Lap Tray, 15" W \$200
Part: MEC1905328-15B. HCPCS: E0950
For use with 16"-17" back width.

- Black Lap Tray, 17" W \$200
Part: MEC1905328-17B. HCPCS: E0950
For use with 18"-19" back width.
- Black Lap Tray, 19" W \$200
Part: MEC1905328-19B. HCPCS: E0950
For use with 20"-22" back width.

1b. Clear Lap Tray

- Clear Lap Tray, 11" W \$200
Part: MEC1905328-11C. HCPCS: E0950
For use with 12"-13" back width.
- Clear Lap Tray, 13" W \$200
Part: MEC1905328-13C. HCPCS: E0950
For use with 14"-15" back width.
- Clear Lap Tray, 15" W \$200
Part: MEC1905328-15C. HCPCS: E0950
For use with 16"-17" back width.

- Clear Lap Tray, 17" W \$200
Part: MEC1905328-17C. HCPCS: E0950
For use with 18"-19" back width.
- Clear Lap Tray, 19" W \$200
Part: MEC1905328-19C. HCPCS: E0950
For use with 20"-22" back width.

2. Mounting Hardware

- Mounting Bracket, Chest Bar/Lap Tray \$75
Part: MEC2104886

23.2. Simple Solutions Lap Tray

Mounting bracket is required with lap tray selection. See chart below for compatibility.

1. Lap Tray

- Simple Solutions Small Tray with Cut Out \$309
Part: ST-SSLXS010001171-Q. HCPCS: E0950
- Simple Solutions Medium Tray with Cut Out \$309
Part: ST-SSLXM010001171-Q. HCPCS: E0950
- Simple Solutions Large Tray with Cut Out \$309
Part: ST-SSLXL010001171-Q. HCPCS: E0950

2. Mounting Hardware

- Swing Away Mounting Hardware \$302
Part: ST-SSLXHD10000060-Q

Simple Solutions Lap Tray Armrest Compatibility	Width											
	10	11	12	13	14	16	17	18	19	20	21	22
Flip-Back	S	S	S	S	S	S	S	M	ML	ML	L	L
Flip-Back, Cane Mounted	S	S	S	S	S	S	M	ML	ML	L	L	
Drop-In, T-Arm	S	S	S	S	S	S	M	ML	ML	L	L	
Two-Post, Flip-Back	S	S	S	S	S	S	M	ML	ML	L	L	
Cantilever	S	S	S	S	S	M	M	ML	ML	L	L	
Q-Captain, Flip-Up						S	M	ML		L		
Single Post/ Heavy Duty/ Bariatric	S	S	S	S	S	S	M	S	ML	ML	L	L

S:ST-SSLXS010001171-Q M:ST-SSLXM010001171-Q L:ST-SSLXL010001171-Q
 :Compatible :Armrest configuration may not hold tray position
 *Some armrests listed may not be available on this order form.

24. ALTERNATIVE CONTROLS AND ELECTRONICS

Optional. Select a type of specialty controls below. More options are available on the [Alternative Controls and Electronics order form](#).

24.1. Enhanced Display

Q-Logic 3 Enhanced Display Kit..... \$1,980
 Part: CTL167726
 Select a color option below. Includes enhanced display, flexible gooseneck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.



Enhanced Display

24.2. Q-Logic SCIM Module

Optional.
SCIM Module for Q-Logic 3 \$525
 Part: CTL2001963
 Allows use of joystick display for specialty devices. Available with Q-Logic 3 EX Joystick selection above. Up to 5 profiles to control drive, seat, and auxiliary functions, excluding sip-n-puff and single switch functions. A switch is required to operate the input device attached to the SCIM module.



SCIM Module

24.3. Q-Logic 3 SAJ

Q-Logic 3 EX SAJ..... \$840
 Part: CTL143980
 Must be used in conjunction with Enhanced Display or Q-Logic 3 Joystick..



EX SAJ

24.4. Head Array

i-Drive 4.0 Bluetooth Tri-array Series, 3 Switch with Egg Switch \$4,774 Part: ST-IDHBT200-U-Q. HCPCS: E2330	Junior, Curved, i-Drive Tri-Array with Mounting Hardware \$4,774 Part: ST-IDHBT450-U-Q. HCPCS: E2330
i-Drive 4.0 Bluetooth Pro Series, 3 Switch with Egg Switch \$4,825 Part: ST-IDHBT300-U-Q. HCPCS: E2330	
Adult, Curved, i-Drive Tri-Array with Mounting Hardware \$4,774 Part: ST-IDHBT400-U-Q. HCPCS: E2330	



Adult Pro-Series

24.5. Sip and Puff Drive

1. Module

Sip and Puff Module for Q-Logic 3..... \$1,990
 Part: CTL168056. HCPCS: E2325
 Includes module and mounting hardware.



Sip & Puff Module

2. Straw and Mount

The straw and mount do not include electrical components to operate the Sip&Puff with a head array. This component requires the selection of a Q-Logic Sip&Puff module or additional module from Stealth Products for operation.

Stealth Sip and Puff for Comfort Series, Left (SP625M-CL) \$318 Part: ST-SP625M-CL-Q. HCPCS: E2326 Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.	Stealth Sip and Puff for Stealth i-Drive, Head Array (SP625-U) \$318 Part: ST-SP625-U-Q. HCPCS: E2326 Includes straw and mounting hardware. Use with Stealth i-Drive Head Arrays.
Stealth Sip and Puff for Comfort Series, Right (SP625M-CR) \$318 Part: ST-SP625M-CR-Q. HCPCS: E2326 Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.	



Stealth Sip & Puff

3. Sip & Puff Accessories

Sip & Puff Straws, 12 pack (SP625-S12).....\$35
 Part: ST-SP625-S12-Q

24.6. Egg Switches

Red Egg Switch (ESRED) \$87 Part: SWTMCRO1077	Green Egg Switch (ESGRN) \$87 Part: SWTMCRO1078
Blue Egg Switch (ESBLU) \$87 Part: SWTMCRO1079	Yellow Egg Switch (ESYEL) \$87 Part: SWTMCRO1080
Black Egg Switch (ESBLK) \$87 Part: SWTMCRO1081	



Egg Switch

Alternative controls and electronics section continued on following page.

24. ALTERNATIVE CONTROLS AND ELECTRONICS CONTINUED

24.7. Switch Adapter

2 to 1 Switch Adapter \$15
Part: DWR1361H037
Connects 2 single switches to a stereo jack to act as a dual switch.

24.8. Environmental Controls

Q-Logic 3 ECU Module Kit \$2,520
Part: CTL168057